# **SPECIFICATIONS**

Particle Size: 2 μm, 2.7 μm

Pore Size: 90 Å

Carbon Load: 6.5% Surface Area:

2 μm: 120 m2/g 2.7 µm: 135 m<sup>2</sup>/g

Endcapped: No Low pH Limit /Max T: 1/90 °C High pH Limit/Max T: 8/40 °C

# PART NUMBERS

2.7 µm ANALYTICAL COLUMNS	
Dimensions: ID x Length (in mm)	Part Number
1.5 x 50	9282X-416
1.5 x 100	9282X-616
1.5 x 150	9282X-716
2.1 x 50	92822-416
2.1 x 100	92822-616
2.1 x 150	92822-716
3.0 x 50	92823-416
3.0 x 100	92823-616
3.0 x 150	92823-716
4.6 x 50	92824-416
4.6 x 100	92824-616
4.6 x 150	92824-716


2.0 μm ANALYTICAL COLUMNS		
Dimensions: ID x Length (in mm)	Part Number	
2.1 x 50	91822-416	
2.1 x 100	91822-616	
2.1 x 150	91822-716	
3.0 x 50	91823-416	
3.0 x 100	91823-616	
3.0 x 150	91823-716	

2.7 μm GUARD COLUMNS		
Guard columns, 3-pack		
Dimensions: ID x Length (in mm)	Part Number	
2.1 x 5	M2822-116	
3.0 x 5	M2823-116	
4.6 x 5	M2824-116	
Guard Column Holder	94900-001	

2.0 µm GUARD COLUMNS  Guard columns, 3-pack		
2.1 x 5	M1822-116	
3.0 x 5	M1823-116	
Guard Column Holder	94900-001	









halocolumns.com

HALO and Fused-Core are registered trademarks of Advanced Materials Technology, Inc.





# HALO® LPH - C18

# **INTRODUCING HALO® LPH-C18**

Introducing a low pH compatible, 90 Å, superficially porous particle C18 phase useful for any chromatographer running under low pH conditions. The sterically protected ligand reduces acidic hydrolysis which enables low pH mobile phases to be used without sacrificing column performance over time.

Si-(CH<sub>2</sub>)<sub>17</sub> — CH<sub>3</sub>

# FEATURES OF HALO® LPH-C18

- Improved stability with low pH mobile phases of pH 1-2
- Highly reproducible alkyl chain bonded phase coverage
- Built upon Fused-Core® Technology for fast, efficient, rugged separations


## **Best Applications:**

Wide range of small molecule applications including:

- polyphenols
- cannabinoids
- pesticides

## QUALITY YOU CAN COUNT ON

A separation of parabens is performed on a HALO 90 Å LPH-C18 column under low pH (pH 1) and high temperature conditions compared to a standard C18 SPP column. Due to the sterically protected ligand, the LPH-C18 column can withstand these conditions and maintain stable retention times.



#### **TEST CONDITIONS**

Column: HALO 90 Å LPH-C18, 2.7 μm 2.1x50 mm

Part Number: 92822-416

Mobile Phase A: Water, 1% TFA (pH: 1)

Mobile Phase B: Acetonitrile
Gradient: Time %B

Time %B 0.0 20 7.50 20 7.51 5 45.00 5 47.00 100 51.00 100

51.01 20

60.00 20

Flow Rate: 0.5 mL/min Pressure: 108 bar Temperature: 60 °C


Detection: UV 254 nm, PDA Injection Volume: 0.4 µL (methyl and ethyl paraben)

Sample Solvent: 25/75 ACN/ Water Data Rate: 100 Hz Response Time: 0.025 sec.

Flow Cell: 1 µl

LC System: Shimadzu Nexera X2

Excellent lot-to-lot reproducible is observed with a mixture of neutral compounds.



#### **TEST CONDITIONS**

Mobile Phase A: Water Mobile Phase B: Acetonitrile Isocratic: 60/40 Acetonitrile/Water

Wavelength: 254 nm

Injection: 2.0 μL (uracil, phenol, 1-Cl-4-nitrobenzene, naphthalene)


Temperature: 30 °C Flow Cell: 1.8 mL/min.

Column: HALO 90 Å LPH-C18 2.7μm 4.6 x 50mm

# **APPLICATIONS**

## COMMON POLYPHENOLS FOUND IN WINE

Common polyphenols found in wine are separated using a HALO 90 Å LPH-C18 column using analytical standards. This stationary phase contains a sterically protected ligand which is ideal for high stability under low



#### **TEST CONDITIONS**

Column: HALO 90 Å LPH-C18, 2.7  $\mu m$  2.1x100 mmMobile Phase A: Water/ 0.1% Formic Acid Mobile Phase B: Acetonitrile/ 0.1% Formic Acid Gradient: Time (min)

0.0	0
3.5	8
7.1	10
25.0	30
26.0	40
27.0	100
29.0	100
30.0	0
35.0	0

Flow Rate: 0.3 mL/min Pressure: 159 bar Temperature: 30 °C Detection: UV 280 nm, PDA Injection Volume: 0.7 µL Sample Solvent: Water Data Rate: 100 Hz Response Time: 0.025 sec.

Flow Cell: 1 µL

LC System: Shimadzu Nexera X2

#### **PEAK IDENTITIES**


- Gallic Acid
- Epigallocatechin 2. 3. Chlorogenic Acid
- 4. Catechin
- 5. Caffeic Acid
- Epicatechin
- Epigallocatechin Gallate
- p-Coumaric Acid Ferulic Acid
- 10. o-Coumaric Acid
- 11. Quercitrin
- 12. Myricetin
- 13. Resveratrol
- Morin

- Quercetin
- 16. Naringenin 17.
- Apigenin

#### 18. Kaempferol


### CATECHINS AND CAFFEINE IN TEA

Catechins belong to the subgroup of polyphenols called flavonoids. These compounds contain antioxidant properties and exist in food and medicinal plants, including tea. An LC-MS separation of catechins and caffeine is demonstrated on a HALO® LPH-C18 column showing excellent resolution using purified standards.



### SEPARATION OF 18 CANNABINOIDS USING HALO® LPH-C18

Anthocyanins, a category of polyphenols, are a type of pigment found in plants that offer several health vegetables, including blueberries. A separation of anthocyanins is performed on a HALO 90 Å LPH-C18 column, which is ideal for the low pH conditions of this method.



#### **TEST CONDITIONS**

Column: HALO 90 Å LPH-C18, 2.7 μm, 4.6 x 150mm Mobile Phase A: 5 mM Ammonium Formate,

0.1% Formic Acid

Mobile Phase B: Acetonitrile, 0.1% Formic Acid

Isocratic: 75% B Flow Rate: 1.5 mL/min Pressure: 232 bar Temperature: 30°C Detection: PDA, UV: 228 nm Injection Volume: 3 uL

Sample Solvent: 75/25 MeOH/ Water

Data Rate: 100 Hz Response Time: 0.025 sec. Flow Cell: 1 µl

LC System: Shimadzu Nexera X2


#### **PEAK IDENTITIES**

- Cannabidivarinic acid (CBDVA)
- Cannabidivarin (CBDV)
- Cannabidiolic acid (CBDA) 3
- Cannabigerolic acid (CBGA)
- Cannabigerol (CBG)
- Cannabidiol (CBD)

- Tetrahydrocannabivarin (THCV)
- Tetrahydrocannabivarinic acid (THCVA)
- Cannabinol (CBN)
- Cannabinolic acid (CBNA)
- Exo-tetrahydrocannabinol (EXO-THC)
- 12. delta 9- Tetrahydrocannabinol (D9-THC)
- 13. delta 8- Tetrahydrocannabinol (D8-THC)
- 14. Cannabicycol (CBL)
- 15. Cannabichromene (CBC)
- 16. Tetrahydrocannabinolic acid A (THCA-A)
- 17. Cannabichromenic acid (CBCA)
- 18. Cannabicyclolic acid (CBLA)

#### BARLEY PESTICIDE SCREENING

Pesticide screening methods can help show whether there is a concern with your soil, crops, and even water supply. A pesticide screening is performed on a sample of barley using a HALO 90 Å LPH-C18 column.



#### **TEST CONDITIONS**

Column: HALO 90 Å LPH-C18 2 µm, 2.1x100 mm Mobile Phase A: Water, 0.1% Formic Acid Mobile Phase B: Acetonitrile, 0.1% Formic Acid

Gradient: Time %B

0.0 30 1.0 30 12.0 100 16.0 100

Flow Rate: 0.2 mL/min Pressure: 235 bar Temperature: 30 °C Detection: +ESI MS/MS Injection Volume: 2 µL Sample Solvent: Methanol

#### **PEAK IDENTITIES**

- 1 Carbendazim
- 2 Dicrotophos
- 3. Azamethiphos
- 4 Pyrimethani 5. Carbofuran
- Dodemorph 6
- Atrazine
- 8. Diuron
- Iprovalicarb 9 Azoxystrobin
- 11 Fluopram
- 12 Methoxyfenozide
- Flutolanil 13.
- 14. Picoxystrobin

