

Chromatography Solutions

Avantor[®] ACE[®], method development kits



Avantor manufactures a range of cutting edge U/HPLC chromatography products from its ISO 9001/ISO 14001 production facility.

The Avantor® ACE® portfolio provides a premium quality product with unique phases that separates what other columns cannot. The range includes novel and traditional stationary phases based on ultra inert silica for excellent reproducibility.



Avantor[®] ACE[®] method development kits

INTELLIGENT SOLUTIONS FOR METHOD DEVELOPMENT

- 5 different ACE Method Development Kits available in a wide range of dimensions for rapid, systematic method development
- Each kit contains carefully selected ACE phases which enables the power of selectivity to be fully exploited
- Each ACE phase provides different selectivity due to differing interactions

FREE

METHOD DEVELOPMENT SUPPORT!

- Not sure which ACE phase or kit will work best for your application?
- FREE Application Support and FREE Method Support Service
- Trust your method development to our experts and free up time for your other projects!

Contact our expert method development team via chromsupport@avantorsciences.com

			SEPARATION MEC	HANISM AND RELA	ATIVE STRENGTH	
	Bonded Phase	Hydrophobic Binding	π-π Interaction	Dipole-Dipole	Hydrogen Bonding	Shape Selectivity
	ACE C18	****	-	-	*	**
ACE Advanced Method Development Kit	ACE C18-AR	****	*** (donor)	*	**	***
(see page 5)	ACE C18-PFP	****	*** (acceptor)	****	***	****
ACE Extended Method Development Kit	ACE SuperC18	****	-	-	-	**
	ACE C18-Amide	****	-	**	****	**/***
see page 9)	ACE CN-ES	***	*	** **** **/** *** **	*	
ACE UltraCore Method	ACE UltraCore SuperC18	***	-	_	-	**
Development Kit see page 13)	ACE UltraCore SuperPhenylHexyl	**	*** (donor)	*	**	***
	ACE C18-300	**	-	-	*	*
ACE Bioanalytical 300Å Method Development Kit	ACE C4-300	*	-	-	-	-
see page 16)	ACE Phenyl-300	*	** (donor)	*	**	**

¹ Approximate value – determined by semi-quantitative mechanism weightings and/or by reference to other ACE phases using >100 characterising analytes.

		SEPARATION MECHANISM AND RELATIVE STRENGTH ²							
			Anionic Analy	te Interactions		Analyte ctions			
	Bonded Phase	Partitioning	Attraction	Repulsion	Attraction	Repulsion	H-bonding		
	ACE HILIC-A	**	-	***	****	-	*		
ACE HILIC Method Development Kit	ACE HILIC-B	***	****	-	-	***	*		
(see page 19)	ACE HILIC-N	****	-	-	-	-	****		

² Approximate value – determined by semi-quantitative mechanism weightings and/or by reference to other ACE phases using >50 characterising analytes.

Why use ACE method development kits?

USING ACE METHOD DEVELOPMENT KITS TO IMPROVE SEPARATIONS

- ACE HPLC/UHPLC columns have earned a well deserved reputation for delivering excellent efficiency, reproducibility and lifetime.
- ACE Method Development Kits group together columns with different mechanisms of interaction to maximise selectivity and improve the likelihood of separating difficult or closely related analytes in mixtures.
- Screening columns containing different bonded phases under the same mobile phase conditions can help you achieve your desired separation more quickly, therefore increasing productivity.

ACE[®] **Stationary Phases** Virtually Eliminate the Negative Effects of Silanols on UHPLC & **HPLC** Separations

COLUMNS WITHIN ACE METHOD DEVELOPMENT KITS PROVIDE ALTERNATIVE SELECTIVITY



ACE C18 provides excellent peak shape, but here the essentially "hydrophobic-only" interaction results in co-elution.

Investigation of alternative bonded phases using the same test conditions is recommended.

ACE C18-PFP provides additional interactions compared to alkyl C18 phases. Whilst these change selectivity, in this instance co-elution of different analyte pairs is observed.

The ACE C18-AR phase provides a further change in selectivity due to different interaction contributions - ultimately enabling a successful separation.

Screening alternative phases can maximise selectivity and reduce method development time.

Sample: 1) metronidazole 2) 4-hydroxybenzoic acid 3) 3-hydroxybenzoic acid 4) benzyl alcohol 5) benzoic acid 6) myrecetin 7) p-cresol 8) propranolol 9) ethyl paraben 10) furosemide 11) anisole 12) 1,3,5-trinitrobenzene 13) toluene 14) nimesulide 15) mefenamic acid 16) 1,2,3-trichlorobenzene Mobile Phase: A = 0.1% formic acid in H2O B = 0.1% formic acid in MeCN - Gradient: 3 - 100% B in 5 minutes Column Dimensions: 50 x 2.1mm - Flow Rate: 0.60ml/min - Temperature: 40°C - Detection: 210nm

ACE Advanced method development kit

- Contains ACE C18, ACE C18-AR and ACE C18-PFP phases

- Ideal starting point for routine method development
- Available in a wide range of dimensions
- Particularly recommended for compounds containing aromatic rings

Phase	Functional Group	Endcapped	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	Recommended pH Range	100% Aqueous Compatible	USP Listing
ACE C18	Octadecyl (C18)	Yes	1.7, 2, 3, 5	100	300	15.5	2.0-8.0°	No	L1
ACE C18-AR	C18 with integral Phenyl	Yes	1.7, 2, 3, 5, 10	100	300	15.5	2.0-8.0°	Yes	L1
ACE C18-PFP	C18 with integral PFP	Yes	1.7, 2, 3, 5, 10	100	300	14.3	2.0-8.0°	Yes	L1

^o For optimum column lifetime, a pH range of 2-8 is recommended. To increase column lifetime at higher pH, organic buffers, low buffer concentrations, high % organic solvent and low temperatures must be considered. Further information is contained within "A Guide to HPLC and LC/MS Buffer Selection" by John Dolan – please contact your distributor to request your FREE copy.

ACE C18	ACE C18-AR	ACE C18-PFP
ACE C18 remains the "go-to" column of choice for HPLC and UHPLC separations. With an excellent reputation for performance, reproducibility and lifetime, ACE C18 provides a rugged, reproducible starting point for method development.	ACE C18-AR combines the excellent performance and advantages of the ACE C18 phase with the added selectivity of an integral phenyl group. Recommended Applications	ACE C18-PFP brings together the stability, reproducibility and low bleed of the ACE C18 phase with the additional selectivity of an integral pentafluorophenyl (PFP) group.
Recommended Applications	 Analytes with π-bonding and conjugated systems 	Recommended Applications - Analytes with π-bonding
 Analytes differing in hydrophobicity Polar, moderately polar and non-polar analytes Uncharged acids and bases Ionized acids or bases using ion-pairing Ideal starting point for method development 	 Analytes with electron delocalization and electron withdrawing groups, such as halogens, nitro groups, ketones, esters and acids Analytes with different dipole moments Analytes differing in hydrophobicity Stereoisomers, steroids, substituted aromatics and sulphur containing compounds 	 Analytes with electron donating groups, such as phenols, aromatic ethers and amines Analytes with proton donor groups Analytes with different dipole moments Analytes differing in hydrophobicity Stuctural isomers, steroids, substituted aromatics and taxanes Fully wettable - 100% aqueous buffer
	 Fully wettable - 100% aqueous buffer compatible Applications where C18 does not provide adequate separation Applications where conventional phenyl phases provide insufficient retention, poor stability, or significant bleed 	 compatible Applications where C18 does not provide adequate separation Applications where conventional PFP phases provide insufficient retention, poor stability or significant bleed

ACE C18 - COMPARISON OF COLUMN INERTNESS

- Column brands from major manufacturers investigated
- Comparison of column efficiency for pyridine a basic molecule



PEAK EFFICIENCY COMPARISON



ACE C18 DELIVERS EXCELLENT PERFORMANCE

RAPID UHPLC SCREENING OF 16 PHARMACEUTICALS AND RELATED COMPOUNDS



Sample: 1) N-acetylprocainamide 2) 3-hydroxybenzoic acid 3) pindolol 4) methylphenylsulphoxide 5) benzyl alcohol 6) quinoxaline 7) 1,4-dinitrobenzene 8) phenacetin 9) 1,2-dimethoxybenzene 10) furosemide 11) anisole 12) methyl benzoate 13) remacemide 14) nimesulide 15) ethyl benzoate 16) diflunisal Mobile Phase: A = 20mM KH₂PO₄, pH 2.7 B = 20mM KH₂PO₄, pH 2.7 in MeOH/H₂O (65:35 v/v) - Gradient: 3 – 100% B in 5 minutes Column Dimensions: 50 x 2.1mm - Flow Rate: 0.60ml/min - Temperature: 60°C - Detection: 214nm.

Application # 1503

LEVERAGING THE UNIQUE SELECTIVITY OF ACE C18-AR

IMPROVING AN ANALGESICS SEPARATION BY CHANGING PHASE



ACE C18-PFP PROVIDES A SEPARATION THAT A C18 OR PFP COLUMN ALONE CANNOT ACHIEVE

THE IMPORTANCE OF MAINTAINING HYDROPHOBICITY DURING MULTI-MODE INTERACTIONS



 Column Dimensions:
 150 x 4.6mm - Sample:
 1) 1,2,3-trimethoxybenzene
 2) 1,2,4-trimethoxybenzene
 3) 1,2-dimethoxybenzene
 4) 1,4-dimethoxybenzene
 5) methoxybenzene

 6) 1,3-dimethoxybenzene
 7) 1,3,5-trimethoxybenzene
 8) neutral molecule (reference)

 Mobile Phase:
 50:50 v/v MeOH/H₂O - Flow Rate:
 1.00ml/min - Temperature:
 40°C - Detection:
 254nm

ACE ADVANCED METHOD DEVELOPMENT KITS

Contains 3 columns: ACE C18, ACE C18-AR and ACE C18-PFP of specified dimensions

	(UHPLC/H	(UHPLC/HPLC hardware format with 1000bar/15000psi pressure limit)					
Column Dimensions	1.7µm	2µm	Зµm	5µm			
2.1 x 50mm	MDKA-17-0502U	MDKA-2-0502U	MDKA-3-0502U	MDKA-5-0502U			
2.1 x 100mm	MDKA-17-1002U	MDKA-2-1002U	MDKA-3-1002U	MDKA-5-1002U			
2.1 x 150mm	-	MDKA-2-1502U	MDKA-3-1502U	MDKA-5-1502U			
2.1 x 250mm	-	-	MDKA-3-2502U	MDKA-5-2502U			
3.0 x 50mm	MDKA-17-0503U	MDKA-2-0503U	MDKA-3-0503U	MDKA-5-0503U			
3.0 x 100mm	MDKA-17-1003U	MDKA-2-1003U	MDKA-3-1003U	MDKA-5-1003U			
3.0 x 150mm	-	MDKA-2-1503U	MDKA-3-1503U	MDKA-5-1503U			
3.0 x 250mm	-	-	MDKA-3-2503U	MDKA-5-2503U			
4.6 x 50mm	-	MDKA-2-0546U	MDKA-3-0546U	MDKA-5-0546U			
4.6 x 100mm	-	MDKA-2-1046U	MDKA-3-1046U	MDKA-5-1046U			
4.6 x 150mm	-	MDKA-2-1546U	MDKA-3-1546U	MDKA-5-1546U			
4.6 x 250mm	-	-	MDKA-3-2546U	MDKA-5-2546U			



ACE Extended method development kit

- Contains ACE SuperC18, ACE C18-Amide and ACE CN-ES phases

- Use ACE SuperC18 to exploit selectivity changes at low, intermediate and high pH
- Available in a wide range of dimensions
- ACE C18-Amide and ACE CN-ES phases both offer alternative selectivity, especially for polar molecules

Phase	Functional Group	Endcapped	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	Recommended pH Range	100% Aqueous Compatible	USP Listing
ACE SuperC18	Octadecyl (C18)	Encapsulated bonding	1.7, 2, 3, 5, 10	90	400	14.8	1.5-11.5°	No	L1
ACE C18-Amide	C18 with integral amide polar group	Yes	1.7, 2, 3, 5, 10	100	300	16.4	2.0-8.0 ^b	Yes	L1/L60
ACE CN-ES	CN with proprietary extended alkyl spacer	Yes	1.7, 2, 3, 5, 10	100	300	12.6	2.0-8.0 ^b	Yes	L10

• ACE SuperC18 is designed for use with LC/MS compatible buffers. Further information is contained within "ACE SuperC18 - A Guide to Buffer Selection" – please contact your distributor to request your FREE copy.

^b For optimum column lifetime, a pH range of 2-8 is recommended. To increase column lifetime at higher pH, organic buffers, low buffer concentrations, high % organic solvent and low temperatures must be considered. Further information is contained within "A Guide to HPLC and LC/MS Buffer Selection" by John Dolan – please contact your distributor to request your FREE copy.

ACE SUPERC18

ACE SuperC18 is a uniquely bonded, EBT endcapped C18 phase which offers unprecedented inertness, excellent efficiency and uncompromising durability over an extended pH range of 1.5 – 11.5.

Recommended Applications

- Analytes differing in hydrophobicityPolar, moderately polar and non-polar
- analytes
- Uncharged acids and bases
- Ionized acids or bases using ion-pairing
- Recommended starting point for developing methods at intermediate and high pH to exploit selectivity changes

ACE C18-AMIDE

ACE C18-Amide is a uniquely designed polar-embedded phase that offers enhanced retention and resolution of polar acidic, phenolic and hydroxysubstituted analytes. The extended spacer ligand technology provides extended column lifetime.

Recommended Applications

- Small water soluble analytes and polar molecules - especially acidic species
- Analytes with H bond donors, acids, bases and phenolic compounds
- Small peptides
- Analytes differing in hydrophobicity
- Fully wettable 100% aqueous buffer compatible
- Applications where C18 does not provide adequate separation
- Applications where conventional amide/polar embedded phases provide insufficient retention, poor stability, or significant bleed

ACE CN-ES

ACE CN-ES is a unique phase having an extended alkyl chain with a terminal cyano group. It provides C18 levels of retention and stability compared to commercial cyano propyl phases which typically exhibit low retentivity and poor stability.

Recommended Applications

- Mixtures of very polar, polar and nonpolar analytes
- Analytes with double and triple bonds
- Analytes differing in hydrophobicity
- Suitable for NP and RP separations
- Fully wettable 100% aqueous buffer compatible
- Applications where a typical C18 column does not provide adequate separation
- Applications where traditional CN bonded phases provide insufficient retention, poor stability or significant bleed
- An orthogonal phase for method development

ACE SUPERC18 PROVIDES EXCELLENT **COLUMN INERTNESS**

- Leading column brands in 50 x 2.1mm LC/MS compatible dimensions at intermediate pH 5.8
- Silica, Hybrid and Superficially Porous particle technologies compared
- Comparison of column efficiency for pyridine a basic molecule
- Efficiency measured at 5% peak height to account for peak tailing effects

PEAK EFFICIENCY COMPARISON

Reproduced with kind permission of The Open University, UK



 $\begin{array}{l} \textbf{Column Dimensions: } 50 \times 2.1 mm (^{\circ} \ 50 \times 2.0 mm) - \textbf{Sample: 1) uracil 2) pyridine 3) phenol - \textbf{Mobile Phase: } 30:70 v/v MeOH/10 mM NH_4OAc in H_2O (pH 5.8) - \textbf{Flow Rate: } 0.20 ml/min - Temperature: 22°C - Detection: 254nm \end{array}$ Comparative data may not be representative of all applications. Application # 1513

Please see back page for acknowledgement of trademark.

USE ACE SUPERC18 TO INVESTIGATE PH EFFECTS

EXPLOIT SELECTIVITY BY ADJUSTING PH



ACE C18-AMIDE PROVIDES ENHANCED POLAR SELECTIVITY

ADVANTAGES OF MULTI-MODE INTERACTIONS FOR HPLC SEPARATIONS

Reproduced with kind permission of The Open University, UK



ACE CN-ES PROVIDES ALTERNATIVE SELECTIVITY

ACE CN-ES PROVIDES A SEPARATION THAT A C18 OR CN COLUMN ALONE CANNOT ACHIEVE



ACE EXTENDED METHOD DEVELOPMENT UHPLC/HPLC COLUMN KITS

Contains 3 columns: ACE SuperC18, ACE C18-Amide and ACE CN-ES of specified dimensions

	(UHPLC/H	PLC hardware format	with 1000bar/15000ps	i pressure limit)
Column Dimensions	1.7µm	2µm	Зµm	5μm
2.1 x 50mm	MDKE-17-0502U	MDKE-2-0502U	MDKE-3-0502U	MDKE-5-0502U
2.1 x 100mm	MDKE-17-1002U	MDKE-2-1002U	MDKE-3-1002U	MDKE-5-1002U
2.1 x 150mm	-	MDKE-2-1502U	MDKE-3-1502U	MDKE-5-1502U
2.1 x 250mm	-	-	MDKE-3-2502U	MDKE-5-2502U
3.0 x 50mm	MDKE-17-0503U	MDKE-2-0503U	MDKE-3-0503U	MDKE-5-0503U
3.0 x 100mm	MDKE-17-1003U	MDKE-2-1003U	MDKE-3-1003U	MDKE-5-1003U
3.0 x 150mm	-	MDKE-2-1503U	MDKE-3-1503U	MDKE-5-1503U
3.0 x 250mm	-	-	MDKE-3-2503U	MDKE-5-2503U
4.6 x 50mm	-	MDKE-2-0546U	MDKE-3-0546U	MDKE-5-0546U
4.6 x 100mm	-	MDKE-2-1046U	MDKE-3-1046U	MDKE-5-1046U
4.6 x 150mm	-	MDKE-2-1546U	MDKE-3-1546U	MDKE-5-1546U
4.6 x 250mm	-	-	MDKE-3-2546U	MDKE-5-2546U



ACE UltraCore method development kit

- Contains ACE UltraCore SuperC18 and SuperPhenylHexyl phases

- Use to exploit selectivity changes at low, intermediate and high pH
- Available in a wide range of dimensions
- Ultra inert core-shell particles and Encapsulated Bonding Technology (EBT^{**}) provide excellent peak shape

Phase	Functional Group	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	Maximum pH Range	USP Listing
ACE UltraCore 2.5 SuperC18	Octadecyl encapsulated	2.5	95	130	7.0	1.5-11.0°	L1
ACE UltraCore 2.5 SuperPhenylHexyl	Phenyl-Hexyl encapsulated	2.5	95	130	4.6	1.5-11.0°	L11
ACE UltraCore 5 SuperC18	Octadecyl encapsulated	5	95	100	5.4	1.5-11.0°	L1
ACE UltraCore 5 SuperPhenylHexyl	Phenyl-Hexyl encapsulated	5	95	100	3.6	1.5-11.0°	L11

• ACE UltraCore columns are designed for use with LC/MS compatible buffers. Further information is contained within "ACE UltraCore – A Guide to Buffer Selection" - please contact your distributor to request your FREE copy.

ACE EXTENDED METHOD DEVELOPMENT UHPLC/HPLC COLUMN KITS

- ACE UltraCore SuperC18 and SuperPhenylHexyl phases are manufactured using our unique Encapsulated Bonding Technology (EBT^{**})
- This technology dramatically increases ligand coverage of the silica surface and effectively eliminates the negative effects of unbonded silanol groups
- The higher ligand coverage results in improved inertness, chromatographic performance and stability

ACE® Stationary Phases Virtually Eliminate the Negative Effects of Silanols on UHPLC & HPLC Separations

ACE ULTRACORE COLUMNS ARE HIGHLY INERT

- Solid-core columns from leading manufacturers investigated
- Comparison of column efficiency for pyridine a basic molecule

PEAK EFFICIENCY COMPARISON

Reproduced with kind permission of The Open University, UK



USE ACE SUPERC18 TO INVESTIGATE PH EFFECTS

EXPLOIT SELECTIVITY BY ADJUSTING PH



INTRODUCING SELECTIVITY CHANGES USING ACE ULTRACORE METHOD DEVELOPMENT KITS

VITAMIN SEPARATION



ACE ULTRACORE METHOD DEVELOPMENT UHPLC/HPLC COLUMN KITS

Contains 2 columns: ACE UltraCore SuperC18 and ACE UltraCore SuperPhenylHexyl of specified dimensions

		ardware format with Opsi pressure limit)
Column Dimensions	2.5µm	5µm
2.1 x 50mm	MDKU-25-0502U	MDKU-5-0502U
2.1 x 100mm	MDKU-25-1002U	MDKU-5-1002U
2.1 x 150mm	MDKU-25-1502U	MDKU-5-1502U
2.1 x 250mm	-	MDKU-5-2502U
3.0 x 50mm	MDKU-25-0503U	MDKU-5-0503U
3.0 x 100mm	MDKU-25-1003U	MDKU-5-1003U
3.0 x 150mm	MDKU-25-1503U	MDKU-5-1503U
3.0 x 250mm	-	MDKU-5-2503U
4.6 x 50mm	MDKU-25-0546U	MDKU-5-0546U
4.6 x 100mm	MDKU-25-1046U	MDKU-5-1046U
4.6 x 150mm	MDKU-25-1546U	MDKU-5-1546U
4.6 x 250mm	-	MDKU-5-2546U



ACE Bioanalytical 300Å method development kit

- Contain ACE C18-300, ACE C4-300 and ACE Phenyl-300 phases

- Ideal starting point for protein and peptide method development
- Available in a wide range of dimensions
- Ultra-inert 300Å phases provide excellent peak shape and reproducibility

Phase	Functional Group	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)	Carbon Load (%)	Recommended pH Range	USP Listing
ACE C18-300	Octadecyl (C18)	3, 5, 10	300	100	9.0	2.0-8.0°	L1
ACE C4-300	Butyl (C4)	3, 5, 10	300	100	2.6	2.0-8.0°	L26
ACE Phenyl-300	Phenyl	3, 5, 10	300	100	5.3	2.0-8.0°	L11

^o For optimum column lifetime, a pH range of 2-8 is recommended. To increase column lifetime at higher pH, organic buffers, low buffer concentrations, high % organic solvent and low temperatures must be considered. Further information is contained within "A Guide to HPLC and LC/MS Buffer Selection" by John Dolan – please contact your distributor to request your FREE copy.

ACE 300Å ULTRA-INERT COLUMNS PROVIDE IMPROVED PEAK SHAPE

ACE 300Å ultra-inert HPLC columns are manufactured using advanced technology that virtually eliminates the negative effects of silanols and metal contamination for the separation of peptides, proteins and other high molecular weight biomolecules. The ultra-inert characteristics of ACE 300Å columns permit the use of as little as 0.005% TFA in the mobile phase. Lower purity columns show unacceptable peak tailing even when using as much as 0.01% TFA. The ability to run at reduced TFA concentrations results in increased sensitivity.

ACE 300Å BIOANALYTICAL COLUMNS PROVIDE EXCELLENT PEAK SHAPE



USING ACE BIOANALYTICAL 300Å METHOD DEVELOPMENT KITS TO OPTIMISE SELECTIVITY

INTRODUCING SELECTIVITY CHANGES BY CAREFUL CONSIDERATION OF BONDED PHASE AND MOBILE PHASE ADDITIVE



The ACE 300Å C18, C4 and Phenyl chemistries contained within ACE Bioanalytical 300Å Method Development Kits enable the analyst to investigate selectivity effects due to phase variations.

The ultra-inert characteristics of the ACE 300Å silica enable different mobile phase additives to be investigated without a deterioration in peak shape or sensitivity.

Columns: 250 x 4.6mm. 5um

Sample: 1) Gly-Tyr 2) Val-Tyr-Val 3) Methionine enkephalin 4) Angiotensin II 5) Leucine enkephalin Mobile Phase: A = 0.1% TFA in H₂O or 0.1% Formic Acid in H₂O (as specified above) B = MeCN Gradient: 10 – 40% B in 25 minutes Flow Rate: 1.00ml/min - Temperature: 22°C - Wavelength: 220nm

Application # 2002

ACE BIOANALYTICAL 300Å METHOD DEVELOPMENT HPLC COLUMN KITS

Contains 3 columns: ACE C18-300, ACE C4-300 and ACE Phenyl-300 of specified dimensions

		lware format with Opsi pressure limit)
Column Dimensions	Зμт	5μm
2.1 x 50mm	MDKB-3-0502	MDKB-5-0502
2.1 x 100mm	MDKB-3-1002	MDKB-5-1002
2.1 x 150mm	MDKB-3-1502	MDKB-5-1502
2.1 x 250mm	-	MDKB-5-2502
3.0 x 50mm	MDKB-3-0503	MDKB-5-0503
3.0 x 100mm	MDKB-3-1003	MDKB-5-1003
3.0 x 150mm	MDKB-3-1503	MDKB-5-1503
3.0 x 250mm	-	MDKB-5-2503
4.6 x 50mm	MDKB-3-0546	MDKB-5-0546
4.6 x 100mm	MDKB-3-1046	MDKB-5-1046
4.6 x 150mm	MDKB-3-1546	MDKB-5-1546
4.6 x 250mm	-	MDKB-5-2546

Guard columns are available for all of these phases

ACE HILIC method development kit

- Contains ACE HILIC-A, ACE HILIC-B and ACE HILIC-N phases
- Alternative and improved selectivity to reversed-phase for polar and very polar analytes
- Available in a wide range of dimensions
- ACE HILIC-A, ACE HILIC-B and ACE HILIC-N provide alternative selectivity to each other

Phase	Functional Group	Endcapped	Particle Size (µm)	Pore Size (Å)	Surface Area (m²/g)		Recommended pH Range	USP Listing
ACE HILIC-A	Proprietary SIL	No	1.7, 3, 5	100	300	-	2.0-7.0	L3
ACE HILIC-B	Proprietary Aminopropyl	No	1.7, 3, 5	100	300	4.0	2.0-7.0	L8
ACE HILIC-N	Proprietary Polyhydroxy	No	1.7, 3, 5	100	300	7.0	2.0-7.0	Pending

WHAT IS HILIC?

- Hydrophilic Interaction Liquid Chromatography (HILIC) was first described by Alpert*
- HILIC is ideal for the separation and retention of polar species including polar neutral and polar ionised analytes
- HILIC separations typically include a polar stationary phase with high organic solvent containing mobile phases
- Mechanistically HILIC is complex (Fig 1) and provides multiple modes of interaction between the analyte, stationary phase, eluent and water enriched layer at the stationary phase particle-eluent interface**

* A. J. Alpert, J. Chromatogr., 499 (1990) 177. ** See the FREE ACE guide to reproducible HILIC method development for more information

WHEN SHOULD YOU CONSIDER HILIC?

- HILIC provides the retention and separation of hydrophilic or polar to very polar analytes not well retained in RPLC
- Hydrophilic or polar to very polar analytes have log P values (measure of lipophilicity) of around zero or less (Fig 2A)
- Generally, polar analytes are suitable for HILIC if they elute before caffeine in gradient RPLC (Fig 2B)



Figure 1. Schematic of interactions between different types of polar analytes and a stationary phase in HILIC mode



Figure 2A. Analyte suitability for HILIC from Log P

Figure 2B. Analyte suitability for HILIC from gradient RPLC

ACE HILIC COLUMNS – 3 ALTERNATIVE SELECTIVITIES

ACE HILIC-A	ACE HILIC-B	ACE HILIC-N
An acidic character phase with an	A basic character phase with an ionisable	A neutral character phase capable of
ionisable negative surface charge	positive surface charge depending on	H-bonding amongst other mechanisms
depending on mobile phase pH	mobile phase pH	of interaction

ACE HILIC-A - An acidic character phase



ACE HILIC-B - A basic character phase





ACE HILIC-N - A polar neutral character phase

ACE HILIC METHOD DEVELOPMENT

- ACE HILIC columns provide alternative selectivity to each other.
- The power of systematic screening of different phase chemistries for HILIC method development is seen below.
- Maximise your HILIC method development success by following the ACE HILIC method development protocol using

three optimised ACE HILIC column chemistries – protocol available in the FREE HILIC Method Development guide.

CONCLUSIONS

ACE HILIC columns provide alternative selectivity to each other – ideal for HILIC method development.

ADVANTAGES OF USING ACE HILIC METHOD DEVELOPMENT KITS



FREE HILIC METHOD DEVELOPMENT TECHNICAL GUIDE

A 38 page HILIC Method Development Technical Guide illustrating a tried and tested approach to HILIC method development is available. Request your copy today and learn how to develop reproducible and robust HILIC methods simply and efficiently.

Alternatively, please contact our technical support team via **chromsupport@avantorsciences.com**

ACE HILIC METHOD DEVELOPMENT UHPLC/HPLC COLUMN KITS

Contains 3 columns: ACE HILIC-A, ACE HILIC-B and ACE HILIC-N of specified dimensions

Column	(UHPLC/HPLC hardware format with 1000bar/15000psi pressure limit)		
Dimensions	1.7µm	Зµm	5µm
2.1 x 50mm	MDKH-17-0502U	MDKH-3-0502U	MDKH-5-0502U
2.1 x 100mm	MDKH-17-1002U	MDKH-3-1002U	MDKH-5-1002U
2.1 x 150mm	-	MDKH-3-1502U	MDKH-5-1502U
2.1 x 250mm	-	MDKH-3-2502U	MDKH-5-2502U
3.0 x 50mm	MDKH-17-0503U	MDKH-3-0503U	MDKH-5-0503U
3.0 x 100mm	MDKH-17-1003U	MDKH-3-1003U	MDKH-5-1003U
3.0 x 150mm	-	MDKH-3-1503U	MDKH-5-1503U
3.0 x 250mm	-	MDKH-3-2503U	MDKH-5-2503U
4.6 x 50mm	-	MDKH-3-0546U	MDKH-5-0546U
4.6 x 100mm	-	MDKH-3-1046U	MDKH-5-1046U
4.6 x 150mm	-	MDKH-3-1546U	MDKH-5-1546U
4.6 x 250mm	-	MDKH-3-2546U	MDKH-5-2546U

UHPLC and HPLC column accessories

UHPLC COLUMN CONNECTORS

- Pressure rating >1700 bar (>25000 psi)
- Compatible with all UHPLC systems¹
- Compatible with all UHPLC column brands
- Eliminates poor connections
- Innovative reusable design

HPLC COLUMN CONNECTORS

- Fingertight to 350 bar (5000 psi)
- Reuseable and simple to install
- Eliminates poor connections
- Compatible with all HPLC column brands and instruments



Reuseable column connector with 10-32 threads for 1/16" od tubing

Proprietary PEEK™ cone

All UHPLC column brands require correct installation in order to realise maximum column efficiency. To avoid connection problems, permanently swaged fittings are not recommended as they do not allow free movement between the tubing, fitting and column inlet on installation. This can result in a poorly connected column that shows unexpected peak tailing due to the introduction of extra column volume (dead volume) to the system. Alternatively, a leak at the inlet fitting connection may be observed.

ACE Excel UHPLC Column Connectors (p/n EXL-CC10, 10 pack) enable the inlet end of UHPLC columns to be correctly installed every time. Their unique reuseable design ensures that they maintain a 1700 bar (25000 psi) pressure rating with repeated use, yet do not permanently swage onto the inlet tubing. To maximise the lifetime of the fitting, the use of a torque wrench (p/n EXL-TW) is required.

At the outlet end of the UHPLC column (where pressure demands are lower but a correct connection remains important), ACE Fingertight HPLC Column Connectors (p/n ACE-CC10, 10 pack, see below) may alternatively be used.

¹Note: For inlet connections onto a Waters Acquity system (containing a Waters Acquity 1/16" fitting and ferrule on the inlet tubing) the use of a pre-column filter incorporating the unique Waters Acquity column port profile is alternatively recommended (p/n EXL-PCF10/ACQ - 10 pack) to ensure maximum compatibility with the Waters Acquity system fittings. ACE Fingertight HPLC Column Connector (p/n ACE-CC10, 10 pack)



ACE Fingertight HPLC Column Connectors (p/n ACE-CC10, 10 pack) are recommended for the connection of both the inlet and outlet ends of HPLC columns.

Manufactured from premium quality PEEK[™], the fittings simply hand tighten to provide a perfect column connection, and are pressure rated to 350 bar/5000 psi.

ACE Fingertight HPLC Column Connectors may additionally be used at the outlet end of UHPLC columns, where pressure demands are lower but a correct connection remains important.



Setting science in motion to create a better world

Avantor[®] is a leading global provider of mission critical products and services to customers in the biopharma, healthcare, education & government, and advanced technologies & applied materials industries. We operate in more than 30 countries and deliver an extensive portfolio of products and services.