


MAbPac SCX-10 Column for Monoclonal Antibody Variant Analysis and Characterization

The Thermo Scientific™ MAbPac™ SCX-10 columns separate closely-related monoclonal antibody variants for characterization and quality control assessment.

- High resolution of monoclonal antibody variants
- Exceptionally high efficiency
- High throughput capability
- Available in 3, 5 or 10 μm particle size
- UHPLC column formats
- Excellent lot-to-lot and column-to-column reproducibility
- No Hydrophobic interactions
- Ideal for R&D and QA/QC development

MAb Variant Resolution Using the MAbPac SCX-10



High-Resolution, High-Efficiency Analysis of Monoclonal Antibody Variants

The MAbPac SCX-10 column is a strong cation-exchange column designed specifically for the high resolution, high efficiency analysis of monoclonal antibodies and associated variants. The unique nonporous pellicular resin provides high resolving power, permitting the separation of monoclonal antibody variants that differ by as little as one charged residue.

The MAbPac SCX-10 columns are available in 10, 5, and 3 μm particle sizes in standard and UHPLC formats. Significantly faster analysis can be achieved on MAbPac SCX-10 columns with smaller particle sizes and in UHPLC formats. The MAbPac SCX-10, 3 and 5 μm particle size columns in a shorter format provide this exceptional capability.

The MAbPac SCX resin technology is the basis for the superior performance of monoclonal antibody variant analysis. The nonporous core particle provides high rates of mass transfer, which results in high efficiency separations. A proprietary hydrophilic layer surrounds the polymeric beads preventing hydrophobic interactions between proteins and the core of the resin. A proprietary grafted cation-exchange surface provides pH selectivity control and high-resolution separations.

The MAbPac SCX-10 column is complementary to the industry-leading Thermo Scientific ProPac WCX-10 column for monoclonal antibody variant analysis, offering an alternative selectivity and providing higher resolution and efficiency for variant analysis of monoclonal antibody samples.

Characterization of Monoclonal Antibody Heterogeneity

Monoclonal antibodies are currently developed by pharmaceutical and biotechnology companies for various therapeutic applications. They undergo several posttranslational modifications including oxidations, deamidations, glycosylation, incomplete C-terminal processing, and others. These modifications cause antibody microheterogeneity or variants. Variations in a monoclonal antibody's composition can impact its activity and stability as a biotherapeutic. Monitoring stability of therapeutic monoclonal antibodies is regarded as essential for demonstrating safety and efficacy of a monoclonal antibody drug, and is expected by the FDA and other regulatory agencies. With its ability to characterize monoclonal heterogeneity, the MAbPac SCX-10 column can be used for stability testing and other monoclonal antibody applications by providing high efficiencies and high resolution of monoclonal antibody variant analysis.

Acidic and Basic Variant Analysis

One of the most important and common analyses of monoclonal antibody heterogeneity is the monitoring and determination of acidic and basic variants. The MAbPac SCX-10 column provides excellent peak efficiencies and high resolution for acidic and basic variant analysis of monoclonal antibodies as shown in Figure 1, where three different monoclonal antibodies (A, B, C) are separated.




Figure 1: The MAbPac SCX-10 column provides excellent peak efficiencies and exceptional resolution of monoclonal antibody variants. A: MAb A separation showing excellent peak efficiencies; B: Separation of MAb B acidic, basic and C-terminal lysine variants; C: MAb C variant separation

C-Terminal Lysine Variant Analysis

During the development and the production of therapeutic monoclonal antibodies, characterization of structural variants is a critical challenge. C-terminal processing of lysine residues on the heavy chain of monoclonal antibodies is a common structural variation that demands analysis. Incomplete monoclonal antibody processing results in charge heterogeneity, which is readily identified using the MAbPac SCX-10 column. Figure 2 illustrates this with the baseline resolution of C-terminal lysine truncation variants, and many other acidic and basic variants of a monoclonal antibody sample. After treatment with carboxypeptidase B, only one major peak remains, verifying that the three major peaks were due to different lysine truncations.

Analysis of MAbs Fragments After Digestion with Papain and Carboxypeptidase

Monitoring stability of therapeutic proteins is regarded as essential for demonstrating safety and efficacy. Sometimes monoclonal antibodies are treated with papain and carboxypeptidase enzymes in order to monitor charge heterogeneity. Monoclonal antibodies treated with papain enzyme are separated into their Fab and Fc fragments. Figure 3 shows the well resolved Fab and Fc fragments after a monoclonal antibody and its variants are treated with papain alone, or with papain and carboxypeptidase together. The expected acidic, C-terminal lysine truncation-containing Fc fragments and Fab fragment peaks are determined and well resolved. Lysine truncation variant peaks collapse into one main peak with the carboxypeptidase treatment.




Figure 2: Baseline resolution of C-terminal lysine variants of a monoclonal antibody sample. A second chromatogram verifies that the three major peaks are due to variations in C-terminal content: after the treatment with carboxypeptidase B, only one major peak remains




Figure 3: High resolution analysis of monoclonal antibody fragments after treatment with papain or papain and carboxypeptidase enzymes. The expected acidic, C-terminal lysine variants containing Fc, and Fab fragment peaks are well resolved

Optimizing the Separation with pH Gradients

High-resolution monoclonal antibody variant separations are achieved through the optimization of different buffers, gradient changes, and pH selectivity control. Monoclonal antibody chromatographic separations can be achieved and optimized using gradients based on changing salt or pH conditions. Figure 4 demonstrates a high-resolution monoclonal antibody variant separation using a pH gradient.

The use of a pH gradient for the analysis of monoclonal antibody samples on MAbPac SCX-10 columns offers some key advantages. A single pH method can be used for the analysis of monoclonal antibody samples with varying iso-electric points. Also, pH gradients can provide high resolution separations on short length 50 mm columns thus providing short separation times or fast analysis. These attributes provide fast high-throughput analysis of multiproduct monoclonal antibody samples.




Figure 4: With the use of a pH gradient, the MAbPac SCX-10 column provides high resolution of monoclonal antibody variants




Figure 5: Fast salt-gradient elution: MAb separations using 3 µm MAb SCX-10 columns

Fast MAB Characterization Analysis

Fast high-resolution analysis of monoclonal antibody samples with exceptionally short separation run times can be achieved on MAbPac SCX-10 columns. This can be accomplished through the use of columns with 3 or 5 μm particle sizes and a column length of 50 mm. Using a 3 μm MAbPac SCX-10 column Figure 5 shows high resolution monoclonal antibody variant separations with exceptionally short separation times which were accomplished with pH and salt gradients, respectively. Figure 6 demonstrates the 5 μm MAbPac SCX-10 column providing fast, high resolution monoclonal antibody variant analysis using MES-based salt gradients. Panel A shows a 5-minute gradient and Panel B shows a 10-minute gradient. High resolution is achieved with a longer, shallow gradient.

The 3 and 5 μm MAbPac SCX-10 columns provide the same high resolution as the MAbPac SCX-10, 10 μm , longer columns but with significantly faster analysis time. Figure 7 shows the 3 μm MAbPac SCX-10 column providing significantly faster analysis while maintaining the comparable high resolution for monoclonal antibody variant analysis as the MAbPac SCX-10, 10 μm , and Thermo Scientific™ ProPac™ WCX-10, longer columns.




Figure 6: The 5 μm MAbPac SCX-10, 4 \times 50 mm column provides fast, high resolution monoclonal antibody variant analysis




Figure 7: The 3 μm MAbPac SCX-10, 4 \times 50 mm column provides significantly faster analysis while maintaining the similar high resolution as that obtained using the MAbPac SCX-10 and ProPac WCX-10, 4 \times 250 mm columns

UHPLC MAb Separation

Higher resolution MAb charge variant UHPLC separations can be achieved using the small particle MAbPac cation-exchange phase, with higher flow rate compatibility, together with specially developed bio-inert PEEK lined stainless steel column hardware. These columns take advantage of smaller resin size as well as longer column length to maximize the resolution of MAb variant separation, and are suitable for operation up to 10,000 Psi. Figure 8 shows an example of MAb separation comparing MAbPac SCX-10, 10 μm , 4 \times 250 mm and UHPLC compatible MAbPac SCX-10, 5 μm , 4.6 \times 250 mm columns. Both columns resolved acidic and basic variants from the main lysine truncation peaks. The resolution values of lysine truncation peaks are better on 5 μm column when compared to the 10 μm column, as expected. Higher pressure compatibility of the column hardware allows to use high flow rates for faster separation (See Figure 9).

A. Column:	MAbPac SCX, 10 μm	B. Column:	MAbPac SCX, 5 μm
Dimension:	4 \times 250 mm	Dimension:	4.6 \times 250 mm
Sample:	MAb 5 mg/mL	Sample:	MAb 5 mg/mL
Injection Volume:	15 μL	Injection Volume:	20 μL
Eluent A:	20 mM MES pH 5.6 + 60 mM NaCl	Eluent A:	20 mM MES pH 5.6 + 60 mM NaCl
Eluent B:	20 mM MES pH 5.6 + 300 mM NaCl	Eluent B:	20 mM MES pH 5.6 + 300 mM NaCl
Flow Rate:	0.76 mL/min	Flow Rate:	1 mL/min
Gradient:	15-40% B in 30 min	Gradient:	15-40% B in 30 min




Figure 8: Comparison MAb Separation on UHPLC, BioRS MAbPac SCX-10, 5 μm with MAbPac SCX-10, 10 μm column




Figure 9: MAb Separation on UHPLC, BioRS MAbPac SCX-10, 5 μm , 4.6 \times 250 mm at 1.5 mL/min flow rate

Reproducibility

The MAbPac SCX-10 column's reproducible manufacturing process eliminates column and lot variability as a concern in methods development and data analysis. MAbPac columns are manufactured and tested under the strictest specifications resulting in unmatched column and lot reproducibility. Figure 10 shows column reproducibility with virtually no change in performance for four different columns. In addition, Figure 11 demonstrates MAbPac SCX-10 lot reproducibility with virtually no change in performance for three different lots.




Figure 10: Demonstration of column-to-column reproducibility with virtually no change in performance for four different MAbPac SCX-10 columns tested




Figure 11: Excellent lot-to-lot reproducibility of the MAbPac SCX-10 column is shown. Columns from each of the three different lots (A, B, C) were used for MAb analysis

Specifications and Operational Parameters

Column construction	PEEK™ and PEEK lined Stainless Steel BioRS columns
Bead size	3 µm, 5 µm, 10 µm
Substrate	Highly crosslinked DVB media
Pellicular layer	Proprietary hydrophilic
Functional group	Sulfonic
Pressure limit	7000 psi for 3 µm, 3000 psi for 10 µm, 10,000 psi for BioRS column
Temperature range	Ambient to 60 °C
pH range	2–12
Dynamic loading capacity	Up to 100 µg, depending on the monoclonal antibody sample (10 µm, 4 × 250 mm)
Typical buffers	MES, or other Good's buffers, tris, phosphate
Solvents	50% acetonitrile if needed for cleaning
Detergent compatibility	Nonionic, anionic, or zwitterionic detergents

Ordering Information

MABPac SCX-10 Analytical Column	Part Number
MABPac SCX-10, 3 µm, Analytical Column (4 × 50 mm)	077907
MABPac SCX-10, 5 µm, Analytical Column (4 × 50 mm)	078656
MABPac SCX-10, 5 µm, Analytical Column (4 × 150 mm)	085198
MABPac SCX-10, 5 µm, Analytical Column (4 × 250 mm)	078655
MABPac SCX-10, 10 µm, Analytical Column (4 × 150 mm)	075602
MABPac SCX-10, 10 µm, Analytical Column (4 × 250 mm)	074625
MABPac SCX-10, 10 µm, Analytical Column SCX-10HT (4 × 50 mm)	075603
MABPac SCX-10, 10 µm, Analytical Column SCX-10 (2 × 250 mm)	075604

UHPLC, BioRS MABPac SCX-10 Analytical Column	Part Number
MABPac SCX-10 RS, 5 µm, Analytical Column (4.6 × 50 mm)	082674
MABPac SCX-10 RS, 5 µm, Analytical Column (4.6 × 150 mm)	085209
MABPac SCX-10 RS, 5 µm, Analytical Column (4.6 × 250 mm)	082673

Lot Select Column Set	Part Number
Lot Select Column Set—Three columns from one resin lot (4 × 250 mm)	SP6864
Lot Select Column Set—One column from each of three resin lots (4 × 250 mm)	SP6865

MABPac SCX-10 Semipreparative Column	Part Number
MABPac SCX-10, 10 µm, Semipreparative Column (9 × 250 mm)	SP6866

MABPac SCX-10 Guard Column	Part Number
MABPac SCX-10 Guard Column (2 × 50 mm)	075749
MABPac SCX-10 Guard Column (4 × 50 mm)	074631

For more information, please visit our website at thermoscientific.com/chromatography

© 2013 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



USA and Canada +1 800 332 3331
France +33 (0)1 60 92 48 34
Germany +49 (0) 2423 9431 20 or 21
United Kingdom +44 (0)1928 534110
Japan +81 3 5826 1615

China +86 21 68654588 +86 10 84193588
 +86 20 83145199 800 810 5118
India +91 22 6742 9494 +91 27 1766 2352
Australia 1 300 735 292 (free call domestic)
New Zealand 0800 933 966 (free call domestic)
All Other Enquiries +44 (0) 1928 534 050

Technical Support
North America +1 800 332 3331
Outside North America +44 (0) 1928 534 440

Thermo
 SCIENTIFIC

Part of Thermo Fisher Scientific