


Product Overview

- Adjustable height from 57 cm to 87 cm (other heights available with UH1/UH4 options)
- 1 liftable Column
- Trespa™ chemical resistant worksurface W45 x D55 cm
- Weight capacity: 150 kg
- Worksurface can easily be lifted up or down by commuting an electrical switch
- Worldwide electrical compatibility
- Recommended use with liquid chromatograph
- HPLC Ionbench is on caster wheels and can easily be moved
- Convenient access to the top of your HPLC or UHPLC
- Safe access to solvents

Liftable

motorized

worksurface



Control :




Standard switch

- Up and Down



Optional position display with 3 memories (ref : BCHPOS)

- Up and Down
- Memory function (3 in total)
- Display for height indication



Solvent waste container (optional ref SOLV)



Slit for pipe and tubing



IEC power strip (optional ref IECKIT)

- including 6 IEC power cords

Photo example :



Detailed Technical specifications :

- Compatible with 230V/50HZ or 110V/60HZ
- Electronic components are all compliant with RoHS directives (Restriction of Hazardous Substances Directive).

External Dimensions	cm	inch
Depth	55 (worksurface depth) 65 (for the lower part of the bench)	21.65 (worksurface depth) 25.59 (for the lower part of the bench)
Width	45 (worksurface width) 55 (for the lower part of the bench)	17.71 (worksurface width) 21.65 (for the lower part of the bench)
Height	From 57 to 87 cm Or with option UH4: From 62 to 97 cm Or with option UH1: From 45 to 75 cm	From 22,4 to 34,3 Or with option UH4 from 24,4 to 37,7 inches Or with option UH1 from 17.7 to 29.5 inches
Weight :	22 kg	48.50 lb

Summary

Description	Part-number
Ionbench-LC BCHLC4555 / 220 Volt	BCHLC4555V220
Ionbench-LC BCHLC4555 / 110 Volt	BCHLC4555V110
Solvent waste container 300x200x200	SOLV
6 IEC plug 220 V 50 Hz, with cord	IECKIT
Optional position display with 3 memories	BCHPOS
Optional Upgraded Height from 24,4 to 37,7 inches (62 to 97 cm)	UH4
Optional Custom Height from 17,7 to 29,5 inches (45 to 75 cm)	UH1