

**LC/MS Foreline Pump Supplies**

Description	Part No.
Oil mist filter kit for E1M18/E2M28	3162-1056
Oil mist cartridge filter for MS40+	G1960-80039
Oil return kit	3162-1057
Foreline pump (rotary pump) oil, Inland 45, 1 L, for E1M18/E2M28	6040-0834
Foreline pump (rotary pump) oil, 4 L, for E1M18/E2M28	6040-0798
Foreline exhaust adapter	59980-20134
Hose clamp	1400-0563
Oil mist filter element for E2M18	1535-4970
KF25 clamp, stainless steel	0100-0549
KF25 coseal (inside clamp)	0100-1597
Exhaust tubing	0890-1727
Pump oil drip pan	G1946-00034



Adding foreline pump (rotary pump) oil, 4 L

**LC/MS Chemicals**

Description	6100 Series Single Quadrupole LC/MS		6300 Series Ion Trap LC/MS		6400 Series Triple Quadrupole LC/MS		6500 Series Accurate-Mass Q-TOF LC/MS		Part No.
ES/APCI positive ion performance standard, 5 x 1 mL ampoules	✓	✓			✓		✓		G2423A
Electrospray LC demo sample (Sulfamix)	✓				✓*				59987-20033
ESI+APCI LC demo sample	✓								G1978-85000
ES negative ion performance standard, 5 x 1 mL ampoules	✓				✓				G2424A
APCI negative ion performance standard, 5 x 1 mL ampoules	✓								G2425A
ES/APCI positive ion performance standard	✓								G1946-85004
Multiple-charge compound performance evaluation sample (horse heart myoglobin)	✓								G2426A
ES/APCI positive ion performance standard, 5 x 1 mL ampoules	✓								G2423A
ES-TOF reference mix, 6 x 2 mL ampoules			✓				✓		G1969-85001
ES-TOF biopolymer reference standard kit			✓						G1969-85003
Flushing solvent	✓	✓	✓	✓	✓	✓	✓		G1969-85026
High purity water, 4 L	✓	✓	✓	✓	✓	✓	✓		8500-2236
Methyl alcohol, 1 L	✓	✓	✓	✓	✓	✓	✓		8500-1867
Ammonium formate	✓	✓	✓	✓	✓	✓	✓		G1946-85021
Formic acid, 5 mL	✓	✓	✓	✓	✓	✓	✓		G2453-85060
Acetonitrile, 1 L	✓	✓	✓	✓	✓	✓	✓		G2453-85050

\*Recommended item for familiarization



031 336 90 00 • www.scantecnordic.se