

VOA Vials - Open Top Closures and Septum

Finneran Products Certified For Science™



- Available in clear or amber borosilicate glass
- Choice of .125" or .100" PTFE/Silicone Septa
- Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with Open Top White Polypropylene Closures and Septum

20mL Vials - Class 1

Cat. No.	Description	Qty
9-107	20mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-107	20mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-106	20mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-106	20mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack

40mL Vials - Class 1

Cat. No.	Description	Qty
9-105	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-105	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-102	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-102	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9-103	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	100 Pack
9A-103	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	100 Pack
9-104	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	144 Pack
9A-104	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	144 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

VOA Vials - Open Top Closures and Septum

Standard

Assembled with Open Top White Polypropylene Closures and Septum

60mL Vials - Class 1

Cat. No.	Description	Qty
9-093	60mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-093	60mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-112	60mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-112	60mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack

Precleaned

Assembled with Open Top White Polypropylene Closures and Septum

20mL Vials - Class 2

Cat. No.	Description	Qty
9-107-2	20mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-107-2	20mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-106-2	20mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-106-2	20mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack

40mL Vials - Class 2

Cat. No.	Description	Qty
9-105-2	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-105-2	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-102-2	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-102-2	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9-103-2	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	100 Pack
9A-103-2	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	100 Pack
9-104-2	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	144 Pack
9A-104-2	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	144 Pack

60mL Vials - Class 2

Cat. No.	Description	Qty
9-093-2	60mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-093-2	60mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-112-2	60mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-112-2	60mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

VOA Vials - Open Top Closures and Septum

Precleaned and Certified

Assembled with Open Top White Polypropylene Closures and Septum

20mL Vials - Class 3

Cat. No.	Description	Qty
9-107-3	20mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-107-3	20mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-106-3	20mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-106-3	20mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack

40mL Vials - Class 3

Cat. No.	Description	Qty
9-105-3	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-105-3	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-102-3	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-102-3	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9-103-3	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	100 Pack
9A-103-3	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	100 Pack
9-104-3	40mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	144 Pack
9A-104-3	40mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	144 Pack

60mL Vials - Class 3

Cat. No.	Description	Qty
9-093-3	60mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9A-093-3	60mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.100" PTFE/Silicone Septa	72 Pack
9-112-3	60mL Clear, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack
9A-112-3	60mL Amber, White Polypropylene 24-414mm Open Top Closure, 0.125" PTFE/Silicone Septa	72 Pack

Low Level TOC Vials - Cleaned and Certified (<10 ppb)

Cat. No.	Description	Qty
9-105-TOC	40mL TOC Certified Clear Vial, 24-414mm Open Top Gray Polypropylene Closure, 0.125" PTFE/Silicone Septa	72 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

VOA Vials - Open Top Closures with Patented Top Hat™ Septa

Finneran Products Certified For Science™



- Available in clear or amber borosilicate glass
- Assembled with the patented Top Hat™ Liner system
- Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with Open Top Black Polypropylene Closures and Top Hat™ Septa

20mL Vials - Class 1

Cat. No.	Description	Qty
9-131	20mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-131	20mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

40mL Vials - Class 1

Cat. No.	Description	Qty
9-130	40mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-130	40mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

60mL Vials - Class 1

Cat. No.	Description	Qty
9-135	60mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-135	60mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

VOA Vials - Open Top Closures with Patented Top Hat™ Septa

Precleaned

Assembled with Open Top Black Polypropylene Closures and Top Hat™ Septa

20mL Vials - Class 2

Cat. No.	Description	Qty
9-131-2	20mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-131-2	20mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

40mL Vials - Class 2

Cat. No.	Description	Qty
9-130-2	40mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-130-2	40mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

60mL Vials - Class 2

Cat. No.	Description	Qty
9-135-2	60mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-135-2	60mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

Precleaned and Certified

Assembled with Open Top Black Polypropylene Closures and Top Hat™ Septa

20mL Vials - Class 3

Cat. No.	Description	Qty
9-131-3	20mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-131-3	20mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

40mL Vials - Class 3

Cat. No.	Description	Qty
9-130-3	40mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-130-3	40mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

60mL Vials - Class 3

Cat. No.	Description	Qty
9-135-3	60mL Clear, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack
9A-135-3	60mL Amber, Black Polypropylene 24-400mm Open Top Closure, 0.125" PTFE/Silicone Top Hat™ Septa	72 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

VOA Vials - Solid Top Closures

Finneran Products Certified For Science™



- Available in clear or amber borosilicate glass
- Assembled with Solid Top PTFE-lined closures
- Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with Solid Top White Polypropylene Closures, PTFE Lined

20mL Vials - Class 1

Cat. No.	Description	Qty
9-121	20mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-121	20mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-088	20mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-088	20mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

40mL Vials - Class 1

Cat. No.	Description	Qty
9-120	40mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-120	40mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-089	40mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-089	40mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

60mL Vials - Class 1

Cat. No.	Description	Qty
9-129	60mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-129	60mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-090	60mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-090	60mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

J.G. Finneran Associates, Inc. — Visit www.JGFinneran.com • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2015 certified guidelines

VOA Vials - Solid Top Closures

Precleaned

Assembled with Solid Top White Polypropylene Closures, PTFE Lined

20mL Vials - Class 2

Cat. No.	Description	Qty
9-121-2	20mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-121-2	20mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-088-2	20mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-088-2	20mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

40mL Vials - Class 2

Cat. No.	Description	Qty
9-120-2	40mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-120-2	40mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-089-2	40mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-089-2	40mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

60mL Vials - Class 2

Cat. No.	Description	Qty
9-129-2	60mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-129-2	60mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-090-2	60mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-090-2	60mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

VOA Vials - Solid Top Closures

Precleaned and Certified

Assembled with Solid Top White Polypropylene Closures, PTFE Lined

20mL Vials - Class 3

Cat. No.	Description	Qty
9-121-3	20mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-121-3	20mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-088-3	20mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-088-3	20mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

40mL Vials - Class 3

Cat. No.	Description	Qty
9-120-3	40mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-120-3	40mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-089-3	40mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-089-3	40mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

60mL Vials - Class 3

Cat. No.	Description	Qty
9-129-3	60mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9A-129-3	60mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	72 Pack
9-090-3	60mL Clear, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack
9A-090-3	60mL Amber, White Polypropylene 24-400mm Solid Top Closure, PTFE Lined	100 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Environmental Vials, Closures and Septum

Finneran Products Certified For Science™



- Available in clear or amber borosilicate glass
- Choice of closure and liner, including patented Top Hat™ liner system
- Fixed (bonded) closures eliminate liner fallout

Vials - Class 1

Cat. No.	Description	Qty
D0376-20	20mL Clear Vial, 28x56mm, 24-400mm Thread	100 Pack
D0390-20	20mL Amber Vial, 28x56mm, 24-400mm Thread	100 Pack
320024-2856	20mL Clear Vial, 28x56mm, 24-400mm Thread	144 Pack
320024-2856A	20mL Amber Vial, 28x56mm, 24-400mm Thread	144 Pack
D0390-30	30mL Amber Vial, 28x72mm, 24-400mm Thread	100 Pack
330024-2872	30mL Clear Vial, 28x72mm, 24-400mm Thread	144 Pack
330024-2872A	30mL Amber Vial, 28x72mm, 24-400mm Thread	144 Pack
D0376-40	40mL Clear Vial, 28x95mm, 24-400mm Thread	100 Pack
D0390-40	40mL Amber Vial, 28x95mm, 24-400mm Thread	100 Pack
340024-2895	40mL Clear Vial, 28x95mm, 24-400mm Thread	144 Pack
340024-2895A	40mL Amber Vial, 28x95mm, 24-400mm Thread	144 Pack
D0376-60	60mL Clear Vial, 27x140mm, 24-400mm Thread	100 Pack
D0390-60	60mL Amber Vial, 27x140mm, 24-400mm Thread	100 Pack
36024-27140	60mL Clear Vial, 27x140mm, 24-400mm Thread	144 Pack
36024-27140A	60mL Amber Vial, 27x140mm, 24-400mm Thread	144 Pack

Closures - Class 1

Cat. No.	Description	Pack	Case
5310-24*	24-400mm Black, Polypropylene Open Hole Closure	100	1,000
5310-24414**	24-414mm Black, Polypropylene Open Hole Closure	100	1,000
5320-24	24-400mm Solid Top, Black Polypropylene Unlined Closure	100	1,000
5320-24W	24-400mm Solid Top, White Polypropylene Unlined Closure	100	1,000
5320-24414	24-414mm Solid Top, Black Polypropylene Unlined Closure	100	1,000
5320-24414W	24-414mm Solid Top, White Polypropylene Unlined Closure	100	1,000
5340-24	24-400mm Black, Polypropylene Glass-Filled Open Hole Closure	100	1,000
5360-24W	24-400mm Solid Top, White Polypropylene, PTFE/F217 Lined	100	1,000

* Available in additional colors; (G) Green, and (W) White

** Available in additional colors; (B) Blue, (G) Green, (GY) Gray, (O) Orange, (PK) Pink, (R) Red, (W) White, and (Y) Yellow

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Environmental Vials, Closures and Septum

Closures with Liners - Class 1

Cat. No.	Description	Qty
38-375/100	24-414mm White, Polypropylene Open Hole Cap, Bonded .060" Thick PTFE/Silicone Liner	100 Pack
34-527	24-414mm, White Polypropylene Open Hole Closure, Fixed 0.125" PTFE/Silicone Liner	100 Pack
34-527/BLK	24-414mm, Black Polypropylene Open Hole Closure, Fixed 0.125" PTFE/Silicone Liner	100 Pack
34-527/144	24-414mm, White Polypropylene Open Hole Closure, Fixed 0.125" PTFE/Silicone Liner	144 Pack
34-527/M	24-414mm, White Polypropylene Open Hole Closure, Fixed 0.125" PTFE/Silicone Liner	1,000 Pack
34-375BLK	24-414mm, Black Polypropylene Open Hole Closure, Fixed 0.060" PTFE/Silicone Liner	100 Pack
D0396-24	24-400mm Black Phenolic Solid Top Closure, PE Cone	100 Pack
D0399-24	24-400mm, White Polypropylene Solid Top Closure, PTFE Lined	100 Pack
D0399-24/144	24-400mm, White Polypropylene Solid Top Closure, PTFE Lined	144 Pack
D0400-24	24-400mm, White Polypropylene Solid Top Closure, F217 Lined	100 Pack
D0400-24/144	24-400mm, White Polypropylene Solid Top Closure, F217 Lined	144 Pack
806550TH-24	24-400mm Black Polypropylene Open Top Closure, .0125" Top Hat™ PTFE/Silicone Lined (Patented)	100 Pack

Septum - Class 1

Cat. No.	Description	Qty
38-140	24mm x 0.015" PTFE	100 Pack
610050-24	24mm x 0.100" PTFE/Silicone Septa, 100 per pack	1,000 Case
612550-24	24mm x 0.125" PTFE/Silicone Septa, 100 per pack	1,000 Case

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Amber Glass Wide Mouth Packers

Finneran Products Certified For Science™



-Amber glass bottles ideal for UV light-sensitive liquid samples

-Choice of closure and liner

-Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Prcleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with White Polypropylene Closure, .015" PTFE Lined

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
9-219	1/2 oz, 15cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-220	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-221	2 oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-222	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-223	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-224	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-225	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-226	80 oz, 2500cc Packer, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

Amber Glass Wide Mouth Packers

Precleaned

Assembled with White Polypropylene Closure, .015" PTFE Lined

Amber Glass Wide Mouth Packers - Class 2

Cat. No.	Description	Qty
9-219-2	1/2 oz, 15cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-220-2	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-221-2	2oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-222-2	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-223-2	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-224-2	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-225-2	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-226-2	80 oz, 2500cc Packer, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

Precleaned & Certified

Assembled with White Polypropylene Closure, .015" PTFE Lined

Amber Glass Wide Mouth Packers - Class 3

Cat. No.	Description	Qty
9-220-3	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-221-3	2 oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-222-3	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-223-3	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-224-3	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-225-3	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Amber Glass Wide Mouth Packers with Assorted Closures

Standard

Assembled with Black Phenolic Closure/PE Cone Lined

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
D0173-1/2	1/2 oz, 15cc Packer, 28-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0173-1	1 oz, 30cc Packer, 28-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0173-2	2 oz, 60cc Packer, 33-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0173-4	4 oz, 120cc Packer, 38-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack

Assembled with White Polypropylene Closure, F217 Lined

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
D0180-1/2	1/2 oz, 15cc Packer, 28-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-1	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-2	2 oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-4	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-8	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0180-16	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0180-32	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0180-80	80 oz, 2500cc Packer, 70-400mm Thread, White Polypropylene Closure, F217 Lined	4 Pack

Assembled White Polypropylene Closure, Unlined

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
D0178-1/2	1/2 oz, 15cc Packer, 28-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-1	1 oz, 30cc Packer, 28-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-2	2 oz, 60cc Packer, 33-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-4	4 oz, 120cc Packer, 38-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-8	8 oz, 250cc Packer, 45-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0178-16	16 oz, 500cc Packer, 53-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0178-32	32 oz, 950cc Packer, 53-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0178-80	80 oz, 2500cc Packer, 70-400mm Thread, White Polypropylene Closure, Unlined	4 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Amber Glass Wide Mouth Packers with Assorted Closures

Standard

Without Closures

Amber Glass Wide Mouth Packers - Class 1		
Cat. No.	Description	Qty
D0179-1/2	1/2 oz, 15cc Packer, 28-400mm Thread	24 Pack
D0179B-1/2	1/2 oz, 15cc Packer, 28-400mm Thread - Bulk Pack	624 Pieces
D0179-1	1 oz, 30cc Packer, 28-400mm Thread	24 Pack
D0179B-1	1 oz, 30cc Packer, 28-400mm Thread - Bulk Pack	432 Pieces
D0179-2	2 oz, 60cc Packer, 33-400mm Thread	24 Pack
D0179-4	4 oz, 120cc Packer, 38-400mm Thread	24 Pack
D0179-8	8 oz, 250cc Packer, 45-400mm Thread	24 Pack
D0179-16	16 oz, 500cc Packer, 53-400mm Thread	12 Pack
D0179-32	32 oz, 950cc Packer, 53-400mm Thread	12 Pack
D0179-80	80 oz, 2500cc Packer, 70-400mm Thread	4 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

J.G. Finneran Associates, Inc. — Visit www.JGFinneran.com • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2015 certified guidelines

Boston Round Glass Bottles

Finneran Products Certified For Science™



- Available in clear or amber borosilicate glass

- Choice of closure and liner

- Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
9-170	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-171	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-172	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9-173	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9-174	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9-175	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
9A-170	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-171	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-172	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-173	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-174	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-175	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

Boston Round Glass Bottles

Precleaned

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Boston Round Bottles - Class 2		
Cat. No.	Description	Qty
9-170-2	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-171-2	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-172-2	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9-173-2	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9-174-2	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9-175-2	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

Assembled with White Polypropylene Closure, .015" PTFE Lined

Amber Glass Boston Round Bottles - Class 2		
Cat. No.	Description	Qty
9A-170-2	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-171-2	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-172-2	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-173-2	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-174-2	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-175-2	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

Precleaned and Certified

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Boston Round Bottles - Class 3		
Cat. No.	Description	Qty
9-170-3	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-171-3	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9-172-3	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9-173-3	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9-174-3	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9-175-3	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Boston Round Glass Bottles

Precleaned and Certified

Assembled with White Polypropylene Closure, .015" PTFE Lined

Amber Glass Boston Round Bottles - Class 3

Cat. No.	Description	Qty
9A-170-3	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-171-3	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-172-3	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Closure, PTFE Lined	24 Pack
9A-173-3	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-174-3	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Closure, PTFE Lined	12 Pack
9A-175-3	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Closure, PTFE Lined	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Boston Round Glass Bottles with Assorted Closures

Standard

Assembled with Black Phenolic Closure/PE Cone Lined

Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0149-1/2	1/2 oz, 15mL Clear Boston Round Bottle, 18-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0149-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0149-2	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0149-4	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0149-8	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack
D0149-16	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack
D0149-32	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack

Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0161-1/2	1/2 oz, 15mL Amber Boston Round Bottle, 18-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0161-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0161-2	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0161-4	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, Black Phenolic Closure, PE Cone Lined	24 Pack
D0161-8	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack
D0161-16	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack
D0161-32	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, Black Phenolic Closure, PE Cone Lined	12 Pack

Assembled with White Polypropylene Closure, F217 Lined

Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0156-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0156-2	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0156-4	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0156-8	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0156-16	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0156-32	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Boston Round Glass Bottles with Assorted Closures

Standard

Assembled with White Polypropylene Closure, F217 Lined

Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0168-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0168-2	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0168-4	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0168-8	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0168-16	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0168-32	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack

Assembled with White Polypropylene Closure, Unlined

Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0154-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0154-2	2 oz, 60mL Clear Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0154-4	4 oz, 125mL Clear Boston Round Bottle, 22-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0154-8	8 oz, 250mL Clear Boston Round Bottle, 24-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0154-16	16 oz, 500mL Clear Boston Round Bottle, 28-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0154-32	32 oz, 1000mL Clear Boston Round Bottle, 33-400mm Thread, White Polypropylene Closure, Unlined	12 Pack

Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0166-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0166-2	2 oz, 60mL Amber Boston Round Bottle, 20-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0166-4	4 oz, 125mL Amber Boston Round Bottle, 22-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0166-8	8 oz, 250mL Amber Boston Round Bottle, 24-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0166-16	16 oz, 500mL Amber Boston Round Bottle, 28-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0166-32	32 oz, 1000mL Amber Boston Round Bottle, 33-400mm Thread, White Polypropylene Closure, Unlined	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Boston Round Glass Bottles with Assorted Closures

Standard

Without Closures

Clear Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0155-1/2	1/2 oz, 15mL Clear Boston Round Bottle, 18-400mm Thread	24 Pack
D0155B-1/2	1/2 oz, 15mL Clear Boston Round Bottle, 18-400mm Thread - Bulk Pack	720 Pieces
D0155-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread	24 Pack
D0155B-1	1 oz, 30mL Clear Boston Round Bottle, 20-400mm Thread - Bulk Pack	432 Pieces
D0155-2	2 oz, 60mL Clear Boston Round Bottle, Cap, 20-400mm Thread	24 Pack
D0155B-2	2 oz, 60mL Clear Boston Round Bottle, Cap, 20-400mm Thread - Bulk Pack	288 Pieces
D0155-4	4 oz, 125mL Clear Boston Round Bottle, Cap, 22-400mm Thread	24 Pack
D0155B-4	4 oz, 125mL Clear Boston Round Bottle, Cap, 22-400mm Thread - Bulk Pack	160 Pieces
D0155-8	8 oz, 250mL Clear Boston Round Bottle, Cap, 24-400mm Thread	12 Pack
D0155-16	16 oz, 500mL Clear Boston Round Bottle, Cap, 28-400mm Thread	12 Pack
D0155-32	32 oz, 1000mL Clear Boston Round Bottle, Cap, 33-400mm Thread	12 Pack

Amber Glass Boston Round Bottles - Class 1

Cat. No.	Description	Qty
D0167-1/2	1/2 oz, 15mL Amber Boston Round Bottle, 18-400mm Thread	24 Pack
D0167B-1/2	1/2 oz, 15mL Amber Boston Round Bottle, 18-400mm Thread - Bulk Pack	720 Pieces
D0167-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread	24 Pack
D0167B-1	1 oz, 30mL Amber Boston Round Bottle, 20-400mm Thread - Bulk Pack	432 Pieces
D0167-2	2 oz, 60mL Amber Boston Round Bottle, Cap, 20-400mm Thread	24 Pack
D0167B-2	2 oz, 60mL Amber Boston Round Bottle, Cap, 20-400mm Thread - Bulk Pack	288 Pieces
D0167-4	4 oz, 125mL Amber Boston Round Bottle, Cap, 22-400mm Thread	24 Pack
D0167B-4	4 oz, 125mL Amber Boston Round Bottle, Cap, 22-400mm Thread- Bulk Pack	160 Pieces
D0167-8	8 oz, 250mL Amber Boston Round Bottle, Cap, 24-400mm Thread	12 Pack
D0167-16	16 oz, 500mL Amber Boston Round Bottle, Cap, 28-400mm Thread	12 Pack
D0167-32	32 oz, 1000mL Amber Boston Round Bottle, Cap, 33-400mm Thread	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Clear Glass Straight Sided Wide Mouth Jars - Short and Tall

Finneran Products Certified For Science™



- Ideal for use in soil, sediment and sludge sampling
- Preassembled with .015" PTFE or F217 lined closed top white polypropylene closures
- Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
9-180	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-181	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-182	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-183	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-184	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
9-190	4 oz, 125mL Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-191	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-192	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-193	64 oz, 2000mL Clear Tall Jar, 83-400mm Thread, White Polypropylene Closure, PTFE Lined	6 Pack
9-194	128 oz, 4000mL Clear Tall Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

Clear Glass Straight Sided Wide Mouth Jars - Short and Tall

Precleaned

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 2

Cat. No.	Description	Qty
9-180-2	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-181-2	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-182-2	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-183-2	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-184-2	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 2

Cat. No.	Description	Qty
9-190-2	4 oz, 125mL Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-191-2	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-192-2	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-193-2	64 oz, 2000mL Clear Tall Jar, 83-400mm Thread, White Polypropylene Closure, PTFE Lined	6 Pack
9-194-2	128 oz, 4000mL Clear Tall Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	4 Pack

Precleaned and Certified

Assembled with White Polypropylene Closure, .015" PTFE Lined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 3

Cat. No.	Description	Qty
9-180-3	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-181-3	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-182-3	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-183-3	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-184-3	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 3

Cat. No.	Description	Qty
9-190-3	4 oz, 125mL Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, PTFE Lined	24 Pack
9-191-3	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-192-3	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, PTFE Lined	12 Pack
9-193-3	64 oz, 2000mL Clear Tall Jar, 83-400mm Thread, White Polypropylene Closure, PTFE Lined	6 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Clear Glass Straight Sided Wide Mouth Jars - Short and Tall with Assorted Closures

Standard

Assembled with White Polypropylene Closure, F217 Lined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0086-2	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0086-4	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0086-8	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0086-16	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0086-32	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0098-1	1 oz, 30mL Clear Tall Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0098-2	2 oz, 60mL Clear Tall Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0098-4	4 oz, 250ml Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
D0098-8	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
D0098-16	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack

Assembled with White Polypropylene Closure, Unlined

Clear Glass Short Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0084-2	2 oz, 60mL Clear Jar, 53-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0084-4	4 oz, 125mL Clear Jar, 58-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0084-8	8 oz, 250mL Clear Jar, 70-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0084-16	16 oz, 500mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0084-32	32 oz, 1000mL Clear Jar, 89-400mm Thread, White Polypropylene Closure, Unlined	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0096-1	1 oz, 30mL Clear Tall Jar, 33-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0096-2	2 oz, 60mL Clear Tall Jar, 38-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0096-4	4 oz, 125mL Clear Tall Jar, 48-400mm Thread, White Polypropylene Closure, Unlined	24 Pack
D0096-8	8 oz, 250mL Clear Tall Jar, 58-400mm Thread, White Polypropylene Closure, Unlined	12 Pack
D0096-16	16 oz, 500mL Clear Tall Jar, 70-400mm Thread, White Polypropylene Closure, Unlined	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Clear Glass Straight Sided Wide Mouth Jars - Short and Tall with Assorted Closures

Standard

Without Closures

Clear Glass Short Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0085-2	2 oz, 60mL Clear Jar, 53-400mm Thread	24 Pack
D0085-4	4 oz, 125mL Clear Jar, 58-400mm Thread	24 Pack
D0085-8	8 oz, 250mL Clear Jar, 70-400mm Thread	24 Pack
D0085-16	16 oz, 500mL Clear Jar, 89-400mm Thread	12 Pack
D0085-32	32 oz, 1000mL Clear Jar, 89-400mm Thread	12 Pack

Clear Glass Tall Straight Sided Wide Mouth Jars - Class 1

Cat. No.	Description	Qty
D0097-1	1 oz, 30mL Clear Tall Jar, 33-400mm Thread	24 Pack
D0097-2	2 oz, 60mL Clear Tall Jar, 38-400mm Thread	24 Pack
D0097-4	4 oz, 125mL Clear Tall Jar, 48-400mm Thread	24 Pack
D0097-8	8 oz, 250mL Clear Tall Jar, 58-400mm Thread	12 Pack
D0097-16	16 oz, 500mL Clear Tall Jar, 70-400mm Thread	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

High Density Polyethylene Wide Mouth Jars

Finneran Products Certified For Science™



HDPE Wide Mouth Jars

- *Preassembled with foamed polyethylene lined polypropylene closures*
- *All HDPE bottles and jars are also available with PTFE-lined closures or unlined closures.*
- *Available in standard, and precleaned in accordance with recommended E.P.A. Protocol*

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Prcleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

Standard

Assembled with White Polypropylene Closure, F217 Lined

HDPE Wide Mouth Jars - Class 1		
Cat. No.	Description	Qty
9-200	2 oz, 60mL Wide Mouth Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-200B	2 oz, 60mL Wide Mouth Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	960 Pieces
9-201	4 oz, 120mL Wide Mouth Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-201B	4 oz, 120mL Wide Mouth Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	540 Pieces
9-202	8 oz, 250mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-202B	8 oz, 250mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	230 Pieces
9-203	16 oz, 500mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
9-203B	16 oz, 500mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	144 Pieces
9-204	32 oz, 1000mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
9-204B	32 oz, 1000mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	72 Pieces
9-205	68 oz, 2000mL Wide Mouth Jar, 100-400mm Thread, White Polypropylene Closure, F217 Lined	6 Pack
9-205B	68 oz, 2000mL Wide Mouth Jar, 100-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	30 Pieces
9-206	128 oz, 4000mL Wide Mouth Jar, 110-400mm Thread, White Polypropylene Closure, F217 Lined	4 Pack
9-206B	128 oz, 4000mL Wide Mouth Jar, 110-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	12 Pieces

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

High Density Polyethylene Wide Mouth Jars

Precleaned

Assembled with White Polypropylene Closure, F217 Lined

HDPE Wide Mouth Jars - Class 2

Cat. No.	Description	Qty
9-200-2	2 oz, 60mL Wide Mouth Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-200B-2	2 oz, 60mL Wide Mouth Jar, 33-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	960 Pieces
9-201-2	4 oz, 120mL Wide Mouth Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-201B-2	4 oz, 120mL Wide Mouth Jar, 38-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	540 Pieces
9-202-2	8 oz, 250mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	24 Pack
9-202B-2	8 oz, 250mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	230 Pieces
9-203-2	16 oz, 500mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
9-203B-2	16 oz, 500mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	144 Pieces
9-204-2	32 oz, 1000mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined	12 Pack
9-204B-2	32 oz, 1000mL Wide Mouth Jar, 53-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	72 Pieces
9-205-2	68 oz, 2000mL Wide Mouth Jar, 100-400mm Thread, White Polypropylene Closure, F217 Lined	6 Pack
9-205B-2	68 oz, 2000mL Wide Mouth Jar, 100-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	30 Pieces
9-206-2	128 oz, 4000mL Wide Mouth Jar, 110-400mm Thread, White Polypropylene Closure, F217 Lined	4 Pack
9-206B-2	128 oz, 4000mL Wide Mouth Jar, 110-400mm Thread, White Polypropylene Closure, F217 Lined - Bulk Pack	12 Pieces

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

J.G. Finneran Associates, Inc. — Visit www.JGFinneran.com • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2015 certified guidelines

High Density Polyethylene Narrow Mouth Bottles

Finneran Products Certified For Science™



HDPE Narrow Mouth Bottles

- Preassembled with foamed polyethylene lined polypropylene closures
- All HDPE bottles and jars are also available with PTFE - lined closures or unlined closures.
- Available in standard, and precleaned in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Pre-cleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

Standard

Assembled with White Polypropylene Closure, F217 Lined

HDPE Narrow Mouth Bottles - Class 1

Cat. No.	Description	Qty
9-210	4 oz, 125mL Narrow Mouth Bottle, 20-410mm Thread, Polypropylene Closure, F217 Lined	24 Pack
9-210B	4 oz, 125mL Narrow Mouth Bottle, 20-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	470 Pieces
9-211	8 oz, 250mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined	24 Pack
9-211B	8 oz, 250mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	230 Pieces
9-212	16 oz, 500mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined	12 Pack
9-212B	16 oz, 500mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	144 Pieces
9-213	32 oz, 1000mL Narrow Mouth Bottle, 28-410mm Thread, Polypropylene Closure, F217 Lined	12 Pack
9-213B	32 oz, 1000mL Narrow Mouth Bottle, 28-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	84 Pieces

Precleaned

Assembled with White Polypropylene Closure, F217 Lined

HDPE Narrow Mouth Bottles - Class 2

Cat. No.	Description	Qty
9-210-2	4 oz, 125mL Narrow Mouth Bottle, 20-410mm Thread, Polypropylene Closure, F217 Lined	24 Pack
9-210B-2	4 oz, 125mL Narrow Mouth Bottle, 20-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	470 Pieces
9-211-2	8 oz, 250mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined	24 Pack
9-211B-2	8 oz, 250mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	230 Pieces
9-212-2	16 oz, 500mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined	12 Pack
9-212B-2	16 oz, 500mL Narrow Mouth Bottle, 24-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	144 Pieces
9-213-2	32 oz, 1000mL Narrow Mouth Bottle, 28-410mm Thread, Polypropylene Closure, F217 Lined	12 Pack
9-213B-2	32 oz, 1000mL Narrow Mouth Bottle, 28-410mm Thread, Polypropylene Closure, F217 Lined - Bulk Pack	84 Pieces

Don't see what you are looking for?
Please contact us. Any combination of products is available.

Laboratory Grade Plastic Bottles

Finneran Products Certified For Science™



- Available in wide mouth or narrow mouth

- Available in natural or amber

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

Standard

Wide Mouth Plastic Laboratory Grade Bottles with Closure

High Density Polyethylene (HDPE), Natural - Class 1

Cat. No.	Description	Qty
D0433-4	125mL Natural Wide Mouth Laboratory Grade Bottle, 38mm Closure	24 Pack
D0433B-4	125mL Natural Wide Mouth Laboratory Grade Bottle, 38mm Closure - Bulk Pack	350 Pack
D0433-8	250mL Natural Wide Mouth Laboratory Grade Bottle, 45mm Closure	24 Pack
D0433B-8	250mL Natural Wide Mouth Laboratory Grade Bottle, 45mm Closure - Bulk Pack	180 Pack
D0433-16	500mL Natural Wide Mouth Laboratory Grade Bottle, 55mm Closure	12 Pack
D0433B-16	500mL Natural Wide Mouth Laboratory Grade Bottle, 55mm Closure - Bulk Pack	100 Pack
D0433-32	1000mL Natural Wide Mouth Laboratory Grade Bottle, 65mm Closure	12 Pack
D0433B-32	1000mL Natural Wide Mouth Laboratory Grade Bottle, 65mm Closure - Bulk Pack	55 Pack

High Density Polyethylene (HDPE), Amber - Class 1

Cat. No.	Description	Qty
D0434-4	125mL Amber Wide Mouth Laboratory Grade Bottle, 38mm Closure	24 Pack
D0434B-4	125mL Amber Wide Mouth Laboratory Grade Bottle, 38mm Closure - Bulk Pack	72 Pack
D0434-8	250mL Amber Wide Mouth Laboratory Grade Bottle, 45mm Closure	24 Pack
D0434B-8	250mL Amber Wide Mouth Laboratory Grade Bottle, 45mm Closure - Bulk Pack	72 Pack
D0434-16	500mL Amber Wide Mouth Laboratory Grade Bottle, 55mm Closure	12 Pack
D0434B-16	500mL Amber Wide Mouth Laboratory Grade Bottle, 55mm Closure - Bulk Pack	48 Pack
D0434-32	1000mL Amber Wide Mouth Laboratory Grade Bottle, 65mm Closure	12 Pack
D0434B-32	1000mL Amber Wide Mouth Laboratory Grade Bottle, 65mm Closure - Bulk Pack	24 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

Laboratory Grade Plastic Bottles

Standard

Narrow Mouth Plastic Laboratory Grade Bottles with Closures

High Density Polyethylene (HDPE), Natural - Class 1

Cat. No.	Description	Qty
D0437-4	125mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure	24 Pack
D0437B-4	125mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	350 Pack
D0437-8	250mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure	24 Pack
D0437B-8	250mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	180 Pack
D0437-16	500mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure	12 Pack
D0437B-16	500mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	100 Pack
D0437-32	1000mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure	12 Pack
D0437B-32	1000mL Natural Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	55 Pack

High Density Polyethylene (HDPE), Amber - Class 1

Cat. No.	Description	Qty
D0438-4	125mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure	24 Pack
D0438B-4	125mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	72 Pack
D0438-8	250mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure	24 Pack
D0438B-8	250mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	72 Pack
D0438-16	500mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure	12 Pack
D0438B-16	500mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	48 Pack
D0438-32	1000mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure	12 Pack
D0438B-32	1000mL Amber Narrow Mouth Laboratory Grade Bottle, 28mm Closure - Bulk Pack	24 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Narrow Mouth Septum Bottles

Finneran Products Certified For Science™



- Preassembled with PTFE/Silicone-lined open top black polypropylene closures

- Available in standard, precleaned, or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Prcleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with Open Top Black Polypropylene (PP) Closures and Septum

Clear Glass Narrow Mouth Septum Bottles - Class 1

Cat. No.	Description	Qty
9-160	1 oz, 30mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-161	2 oz, 60mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-162	4 oz, 125mL Clear Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-163	8 oz, 250mL Clear Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Amber Glass Narrow Mouth Septum Bottles - Class 1

Cat. No.	Description	Qty
9A-160	1 oz, 30mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-161	2 oz, 60mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-162	4 oz, 125mL Amber Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-163	8 oz, 250mL Amber Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

Narrow Mouth Septum Bottles

Precleaned

Assembled with Open Top Black Polypropylene (PP) Closures and Septum

Clear Glass Narrow Mouth Septum Bottles - Class 2

Cat. No.	Description	Qty
9-160-2	1 oz, 30mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-161-2	2 oz, 60mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-162-2	4 oz, 125mL Clear Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-163-2	8 oz, 250mL Clear Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Amber Glass Narrow Mouth Septum Bottles - Class 2

Cat. No.	Description	Qty
9A-160-2	1 oz, 30mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-161-2	2 oz, 60mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-162-2	4 oz, 125mL Amber Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-163-2	8 oz, 250mL Amber Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Precleaned and Certified

Assembled with Open Top Black Polypropylene (PP) Closures and Septum

Clear Glass Narrow Mouth Septum Bottles - Class 3

Cat. No.	Description	Qty
9-160-3	1 oz, 30mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-161-3	2 oz, 60mL Clear Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-162-3	4 oz, 125mL Clear Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9-163-3	8 oz, 250mL Clear Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

Amber Glass Narrow Mouth Septum Bottles - Class 3

Cat. No.	Description	Qty
9A-160-3	1 oz, 30mL Amber Septum Bottle, 20-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-161-3	2 oz, 60mL Amber Septum Bottle, 20-40mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-162-3	4 oz, 125mL Amber Septum Bottle, 22-400mm Thread, Black PP Closure, PTFE/Silicone Liner	24 Pack
9A-163-3	8 oz, 250mL Amber Septum Bottle, 24-400mm Thread, Black PP Closure, PTFE/Silicone Liner	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Wide Mouth Septum Jars

Finneran Products Certified For Science™



-Preassembled with ultrasonically welded PTFE/Silicone liner

-Available in standard, precleaned or precleaned/certified in accordance with recommended E.P.A. Protocol

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

CLASS 2: (Precleaned) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried, and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a copy of wash "Standard Operating Procedure."

CLASS 3: (Certified) Containers are processed according to EPA recommended wash procedures. Washed with laboratory grade biodegradable, non-phosphate detergent, rinsed with ASTM Type I De-ionized water, oven dried and assembled in a contaminate-free environment. The exterior of each case is labeled with a lot number. Each case contains a "Certificate of Analysis."

An independent laboratory is used to certify the container's cleanliness, assuring impartial results.

Standard

Assembled with Open Top White Polypropylene Closures and Septum

Tall Straight Sided Wide Mouth Septum Jars - Class 1

Cat. No.	Description	Qty
9-141	4 oz, 125mL Clear Tall Septum Jar, 48-400mm Closure, PTFE/Silicone Liner	24 Pack
9-142	8 oz, 250mL Clear Tall Septum Jar, 58-400mm Closure, PTFE/Silicone Liner	12 Pack

Short Straight Sided Wide Mouth Septum Jars - Class 1

Cat. No.	Description	Qty
9-140	2 oz, 60mL Clear Short Septum Jar, 53-400mm Closure, PTFE/Silicone Liner	24 Pack

Amber Glass Wide Mouth Packers - Class 1

Cat. No.	Description	Qty
9-150	500cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack
9-151	950cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack

Don't see what you are looking for?
Please contact us. Any combination of products is available.

Wide Mouth Septum Jars

Precleaned

Assembled with Open Top White Polypropylene Closures and Septum

Tall Straight Sided Wide Mouth Septum Jars - Class 2

Cat. No.	Description	Qty
9-141-2	4 oz, 125mL Clear Tall Septum Jar, 48-400mm Closure, PTFE/Silicone Liner	24 Pack
9-142-2	8 oz, 250mL Clear Tall Septum Jar, 58-400mm Closure, PTFE/Silicone Liner	12 Pack

Short Straight Sided Wide Mouth Septum Jars - Class 2

Cat. No.	Description	Qty
9-140-2	2 oz, 60mL Clear Short Septum Jar, 53-400mm Closure, PTFE/Silicone Liner	24 Pack

Amber Glass Wide Mouth Packers - Class 2

Cat. No.	Description	Qty
9-150-2	500cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack
9-151-2	950cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack

Precleaned and Certified

Assembled with Open Top White Polypropylene Closures and Septum

Tall Straight Sided Wide Mouth Septum Jars - Class 3

Cat. No.	Description	Qty
9-141-3	4 oz, 125mL Clear Tall Septum Jar, 48-400mm Closure, PTFE/Silicone Liner	24 Pack
9-142-3	8 oz, 250mL Clear Tall Septum Jar, 58-400mm Closure, PTFE/Silicone Liner	12 Pack

Short Straight Sided Wide Mouth Septum Jars - Class 3

Cat. No.	Description	Qty
9-140-3	2 oz, 60mL Clear Short Septum Jar, 53-400mm Closure, PTFE/Silicone Liner	24 Pack

Amber Glass Wide Mouth Packers - Class 3

Cat. No.	Description	Qty
9-150-3	500cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack
9-151-3	950cc Amber Packer, 53-400mm Closure, PTFE/Silicone Liner	12 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Wide Mouth Straight Sided Polypropylene Plastic Jars

Finneran Products Certified For Science™



Standard

- Closures available with additional lining materials

CLASS 1: (Standard) Containers are assembled with liner and closure without washing treatment.

Assembled with White Polypropylene Closure, F217 Lined

Wide Mouth Straight Sided Polypropylene Jars - Class 1

Cat. No.	Description	Qty
D0044-1	1 oz, 30mL Polypropylene Jar, 43-400mm Closure, F217 Lined	24 Pack
D0044B-1	1 oz, 30mL Polypropylene Jar, 43-400mm Closure, F217 Lined - Bulk Pack	1144 Pieces
D0044-2	2 oz, 60mL Polypropylene Jar, 53-400mm Closure, F217 Lined	24 Pack
D0044B-2	2 oz, 60mL Polypropylene Jar, 53-400mm Closure, F217 Lined - Bulk Pack	420 Pieces
D0044-4	4 oz, 125mL Polypropylene Jar, 58-400mm Closure, F217 Lined	24 Pack
D0044B-4	4 oz, 125mL Polypropylene Jar, 58-400mm Closure, F217 Lined - Bulk Pack	324 Pieces
D0044-8	8 oz, 250mL Polypropylene Jar, 70-400mm Closure, F217 Lined	24 Pack
D0044B-8	8 oz, 250mL Polypropylene Jar, 70-400mm Closure, F217 Lined - Bulk Pack	175 Pieces
D0044-16	16 oz, 500mL Polypropylene Jar, 89-400mm Closure, F217 Lined	12 Pack
D0044B-16	16 oz, 500mL Polypropylene Jar, 89-400mm Closure, F217 Lined - Bulk Pack	140 Pieces
D0044-32	32 oz, 950mL Polypropylene Jar, 120-400mm Closure, F217 Lined	12 Pack
D0044B-32	32 oz, 950mL Polypropylene Jar, 120-400mm Closure, F217 Lined - Bulk Pack	48 Pieces

Assembled with White Polypropylene Closure, Unlined

Wide Mouth Straight Sided Polypropylene Jars - Class 1

Cat. No.	Description	Qty
D0045-1	1 oz, 30mL Polypropylene Jar, 43-400mm Closure, Unlined	24 Pack
D0045B-1	1 oz, 30mL Polypropylene Jar, 43-400mm Closure, Unlined - Bulk Pack	1144 Pieces
D0045-2	2 oz, 60mL Polypropylene Jar, 53-400mm Closure, Unlined	24 Pack
D0045B-2	2 oz, 60mL Polypropylene Jar, 53-400mm Closure, Unlined - Bulk Pack	420 Pieces
D0045-4	4 oz, 125mL Polypropylene Jar, 58-400mm Closure, Unlined	24 Pack
D0045B-4	4 oz, 125mL Polypropylene Jar, 58-400mm Closure, Unlined - Bulk Pack	324 Pieces
D0045-8	8 oz, 250mL Polypropylene Jar, 70-400mm Closure, Unlined	24 Pack
D0045B-8	8 oz, 250mL Polypropylene Jar, 70-400mm Closure, Unlined - Bulk Pack	175 Pieces
D0045-16	16 oz, 500mL Polypropylene Jar, 89-400mm Closure, Unlined	12 Pack
D0045B-16	16 oz, 500mL Polypropylene Jar, 89-400mm Closure, Unlined - Bulk Pack	140 Pieces
D0045-32	32 oz, 950mL Polypropylene Jar, 120-400mm Closure, Unlined	12 Pack
D0045B-32	32 oz, 950mL Polypropylene Jar, 120-400mm Closure, Unlined - Bulk Pack	48 Pieces

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

PTFE Cap Liners and Solid Top Closures

Finneran Products Certified For Science™



- Solid virgin PTFE liners can be manually inserted into closure and maintains a snug friction fit.

- PTFE liners are also available in other sizes and thicknesses

PTFE Cap Liners

PTFE Cap Liners .015" Thick

Cat. No.	Description	Qty
38-134	18mm PTFE Liner	100 Pack
38-136	20mm PTFE Liner	100 Pack
38-138	22mm PTFE Liner	100 Pack
38-140	24mm PTFE Liner	100 Pack
38-142	28mm PTFE Liner	100 Pack
38-144	33mm PTFE Liner	100 Pack
38-146	38mm PTFE Liner	100 Pack
38-148	43mm PTFE Liner	100 Pack
38-150	45mm PTFE Liner	100 Pack
38-152	48mm PTFE Liner	100 Pack
38-154	53mm PTFE Liner	100 Pack
38-156	58mm PTFE Liner	100 Pack
38-158	63mm PTFE Liner	100 Pack
38-160	70mm PTFE Liner	100 Pack
38-162	83mm PTFE Liner	100 Pack
38-164	89mm PTFE Liner	100 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

PTFE Cap Liners and Solid Top Closures



Solid Top Closures

- Closures available with various lining materials

Phenolic Closures - Polycone Lined

Cat. No.	Description	Qty
D0396-13	13-425mm Phenolic Cap/PE Cone Lined	100 Pack
D0396-15	15-425mm Phenolic Cap/PE Cone Lined	100 Pack
D0396-18	18-400mm Phenolic Cap/PE Cone Lined	100 Pack
D0396-20	20-400mm Phenolic Cap/PE Cone Lined	100 Pack
D0396-22	22-400mm Phenolic Cap/PE Cone Lined	100 Pack
D0396-24	24-400mm Phenolic Cap/PE Cone Lined	100 Pack
D0396-28	28-400mm Phenolic Cap/PE Cone Lined	100 Pack
D0396-33	33-400mm Phenolic Cap/PE Cone Lined	100 Pack
D0396-38	38-400mm Phenolic Cap/PE Cone Lined	100 Pack

Polypropylene Closures - PTFE Lined

Cat. No.	Description	Qty
D0399-13	13-425mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-15	15-425mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-18	18-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-20	20-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-22	22-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-24	24-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-28	28-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-33	33-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-38	38-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-43	43-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-45	45-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-48	48-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-53	53-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-58	58-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-63	63-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-70	70-400mm White Polypropylene Closure/PTFE Lined	100 Pack
D0399-89	89-400mm White Polypropylene Closure/PTFE Lined	100 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

PTFE Cap Liners and Solid Top Closures

Solid Top Closures

Polypropylene Closures - F217 Lined

Cat. No.	Description	Qty
D0400-18	18-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-20	20-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-22	22-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-24	24-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-28	28-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-33	33-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-38	38-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-43	43-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-45	45-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-48	48-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-53	53-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-58	58-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-63	63-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-70	70-400mm White Polypropylene Closure/F217 Lined	100 Pack
D0400-89	89-400mm White Polypropylene Closure/F217 Lined	100 Pack

Polypropylene Closures - Unlined

Cat. No.	Description	Qty
D0402-18	18-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-20	20-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-22	22-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-24	24-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-28	28-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-33	33-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-38	38-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-43	43-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-45	45-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-48	48-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-53	53-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-58	58-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-63	63-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-70	70-400mm White Polypropylene Closure/Unlined	100 Pack
D0402-89	89-400mm White Polypropylene Closure/Unlined	100 Pack

*Don't see what you are looking for?
Please contact us. Any combination of products is available.*

Standard EPA Vial Box Configuration

Finneran Products Certified For Science™



72 pack: 20 & 40ml vials



100 pack: 20 & 40ml vials



144 pack: 20, 40 & 60ml vials



72 pack: 60ml vials

Other packaging configurations may be available.
Please contact Customer Service.



Analys & Mätteknik

031 336 90 00 • www.scantecnordic.se

J.G. Finneran Associates, Inc. — Visit www.JGFinneran.com • 1-800-552-3696 • 1-856-696-3605

Engineered and manufactured in compliance with ISO 9001:2015 certified guidelines