

Kromasil Eternity columns

Kromasil Eternity, 2.1 mm i.d. columns

Family	Phase	particle size [μm]	column size, i.d. \times length [mm]		
			2.1 \times 50	2.1 \times 100	2.1 \times 150
Eternity	C18	2.5	EH2CLD05	EH2CLD10	
Eternity	C18	5	E05CLD05		E05CLD15
Eternity	PhenylHexyl	2.5	EH2PXD05	EH2PXD10	
Eternity	PhenylHexyl	5	E05PXD05		E05PXD15
EternityXT	C18	1.8	XF1CLD05	XF1CLD10	
EternityXT	C18	2.5	XH2CLD05	XH2CLD10	
EternityXT	C18	5	X05CLD05		X05CLD15
EternityXT	PhenylHexyl	1.8	XF1PXD05	XF1PXD10	
EternityXT	PhenylHexyl	2.5	XH2PXD05	XH2PXD10	
EternityXT	PhenylHexyl	5	X05PXD05		X05PXD15

Kromasil Eternity, 4.6 mm i.d. columns

Family	Phase	particle size [μm]	column size, i.d. \times length [mm]			
			4.6 \times 50	4.6 \times 100	4.6 \times 150	4.6 \times 250
Eternity	C18	2.5	EH2CLA05	EH2CLA10		
Eternity	C18	5	E05CLA05	E05CLA10	E05CLA15	E05CLA25
Eternity	PhenylHexyl	2.5	EH2PXA05	EH2PXA10		
Eternity	PhenylHexyl	5	E05PXA05	E05PXA10	E05PXA15	E05PXA25
EternityXT	C18	2.5	XH2CLA05	XH2CLA10		
EternityXT	C18	5	X05CLA05	X05CLA10	X05CLA15	X05CLA25
EternityXT	C18	10				X10CLA25
EternityXT	PhenylHexyl	2.5	XH2PXA05	XH2PXA10		
EternityXT	PhenylHexyl	5	X05PXA05	X05PXA10	X05PXA15	X05PXA25

Kromasil Eternity, 10 mm i.d. columns

Family	Phase	particle size [µm]	column size, i.d. × length [mm]		
			10 × 100	10 × 150	10 × 250
Eternity	C18	5	E05CLP10	E05CLP15	E05CLP25
Eternity	PhenylHexyl	5	E05PXP10	E05PXP15	E05PXP25
EternityXT	C18	5	X05CLP10	X05CLP15	X05CLP25
EternityXT	C18	10	X10CLP10	X10CLP15	X10CLP25
EternityXT	PhenylHexyl	5	X05PXP10	X05PXP15	X05PXP25

Kromasil Eternity, 21.2 mm i.d. columns

Family	Phase	particle size [µm]	column size, i.d. × length [mm]		
			21.2 × 100	21.2 × 150	21.2 × 250
Eternity	C18	5	E05CLQ10	E05CLQ15	E05CLQ25
Eternity	PhenylHexyl	5	E05PXQ10	E05PXQ15	E05PXQ25
EternityXT	C18	5	X05CLQ10	X05CLQ15	X05CLQ25
EternityXT	C18	10	X10CLQ10	X10CLQ15	X10CLQ25
EternityXT	PhenylHexyl	5	X05PXQ10	X05PXQ15	X05PXQ25

Kromasil Eternity, 30 mm i.d. columns

Family	Phase	particle size [µm]	column size, i.d. × length [mm]		
			30 × 100	30 × 150	30 × 250
Eternity	C18	5	E05CLR10	E05CLR15	E05CLR25
Eternity	PhenylHexyl	5	E05PXR10	E05PXR15	E05PXR25
EternityXT	C18	10	X10CLR10	X10CLR15	X10CLR25

