

Robust and easy to use, this model offers ergonomic shape and lightweight. It provides smooth control in both aspirating and dispensing, using glass or plastic pipettes of volume up to 100 mL. One-year warranty.

Profiller™ manual 437

Product advantages

- Large capacity bulb
- Soft action lever button
- Efficient blow-out system
- PTFE protection membrane filter
- Maintenance free
- Alternative colour-coded nose cones

Squeezable bulb ①

Simply squeeze large capacity silicone bulb for efficient aspiration. If additional aspiration is needed, it can be squeezed even during pipetting.



Precise control ②

Thumb lever button for left and right-handed users to easily control both aspiration and dispensing speeds.

Liquid blow-out ③

Press large button to expel residual liquid when using blow-out type pipettes.

Ordering information - Profiller™ manual 437

Description	Packaging	Cat. No.
Profiller™ manual 437	1 / pk	437.100

437

Pipettes up to 100 mL



The Profillette™ 406

Product advantages

- One-hand operation
- Easy filling and dispensing
- Air inlet valve for fast distribution
- Colour coding per size
- Maintenance free



406

2 mL
10 mL
25 mL



Ordering information - Profillette™ 406

Volume	Colours	Packaging	Cat. No.
2 mL	blue	1 / pk	406.002
10 mL	green	1 / pk	406.010
25 mL	red	1 / pk	406.025
2 / 10 / 25 mL	assorted	3 / pk	406.300

Protection filter ④

Interchangeable hydrophobic membrane filter guarantees effective instrument and sample protection.

Colour coding cones ⑤

Nose cone can be changed for laboratory or department identification. Choose among regular white or alternative blue and green supplied.

Pipette holder

Conical shape of silicone pipette holder tightly fits graduated glass or plastic pipettes. Autoclavable.

Ordering information - Accessories

Description	Packaging	Cat. No.
Profillette™ 406		
Spare silicone pipette holder	1 / pk	1.406.01
Replacement nose screw	1 / pk	1.406.02
Profiller™ manual 437		
Membrane filter, 0.45 µm	5 / pk	322.447
Silicone pipette holder	1 / pk	1.437.01
Nose cone - white	3 / pk	1.437.02
Silicone filter holder	1 / pk	1.437.03