

Desk for particularly heavy equipment

- An extra-strong solution for storing and transferring heavy equipment.
- The desk's reinforced structures make it extremely sturdy in all directions.
- Special castors make it stable and enable it to be moved, for example during the maintenance of the device it is carrying.
- The simple design makes it easy to clean.
- The under-desk furniture can be fitted onto the desk frame, off the floor. This solution allows furniture to be moved horizontally to any location under the desk. The furniture can also be removed and placed under another desk easily.
- Like all our workbenches, this desk is also available with electrical height adjustment. See more information about electric desks [here](#).
- For Standard solutions we have selected high-quality and most commonly used options in terms of dimensions and materials.
- Custom solutions are always tailored to the customer's needs.
- Note the figures in the table below indicate the worktop dimensions.



Standard

Description	Product code and standard depth [mm]	Width [mm]	Height, mm
VISU equipment table	RHE-W-600-H	500-2100	750/900
	RHE-W-800-H	500-2100	750/900
VISU equipment table, electric height adjustment	RLE-W-600-S	500-2100	Adjustment range from 670 to 1320
	RLE-W-800-S	500-2100	Adjustment range from 670 to 1320

Custom


In addition to Standard solutions, we also implement workstations according to customer's wishes. In Custom solutions, dimensions, materials and equipment are completely up to the customer.

Feel free to contact our experts and we will find suitable solutions together!


Features

Installation	Free-standing on the floor, adjustable pads as standard.
---------------------	--

VISU KALUSTE OY

 Eeliksentie 6, 34870 Visuvesi, Finland

 visu@visu.fi

 +358 44 749 4940


 visu.fi




INNOSTAVIA YMPÄRISTÖJÄ
TIETEELLE JA TUTKIMUKSELLE

Load-bearing capacity	At least 240 kg or as necessary.
Materials	Powder coated steel. Also available in stainless steel.
Worktop	Our collection offers a number of options for worktop material. Read more information under worktop materials .
Castors	The furniture can be equipped with special castors that also act as fully lockable adjustable legs.
By order	The frame is also available with a lower level.

VISU KALUSTE OY

 Eeliksentie 6, 34870 Visuvesi, Finland

 visu@visu.fi

 +358 44 749 4940

 visu.fi